B.Sc. CHEMISTRY LAB MANUAL 3rd Semester

Prepared By Pure & Applied Science Dept. Chemistry

MIDNAPORE CITY COLLEGE

CHEMISTRY HONOURS [Choice Based Credit System] (Semester Programme) SEMESTER-III

C5P: Physical Chemistry Lab

Experiment 1: Study of viscosity of unknown liquid (glycerol, sugar) with respect to water

Theory:

From the literature value of viscosity coefficient of water (η_w) at the experimental temperature and the measured specific gravity of solutions of different concentrations we can determine the value of absolute viscosity coefficients of experimental solution of different concentrations. Thus we can study the variation of viscosity coefficient with concentration by plotting a graph of viscosity coefficient vs. concentration. We can also determine concentration of a solution of unknown concentration from the calibration curve.

Apparatus Required:

- 1) 100ml beaker -2
- 2) 250ml volumetric flask -1
- 3) 50ml volumetric flask -4
- 4) 10ml sp. Gravity bottle -1
- 5) Viscometer (Time of flow of 10ml of water should be at least 80 secs)
- 6) 10ml pipette -1
- 7) Stop watch.

Chemical required:

Sucrose/ Glycerol

Procedure:

- 1) Prepare a stock solution of 15% (v/v) Glycerol in a 250ml volumetric flask. Then prepare 12%, 9%, 6%, 3% Glycerol solution each of 50ml from the stock solution by exact dilution.
- 2) Clean the viscometer with chromic acid and wash thoroughly with deionized water. Remove the water completely. Add 10 ml water in the wider limb using a 10ml pipette. Suck up the water in the other limb and allow it to run between two specified marks. Note the time of flow. Remove the water completely. Rinse the viscometer with experimental liquid (start with lowest concentration of experimental solution and then with increasing concentrations] and discard the rinsing. Add 10ml of experimental liquid with pipette , repeat the procedure and note the time

of flow. Note the time of flow for each of the liquid (water, experimental liquids) at least twice.

3) Use a clean dry 10ml specific gravity bottle to determine the specific gravity of the experimental liquids at room temperature. Record the temperature.

Experimental Result :

- 1) Room temperature:
- 2) Preparation of Glycerol solution of different concentration from the stock (15%) solution.

Concentration of the prepared solution (%)	Volume of 15% Glycerol solution(ml)	Volume of water(ml)
3	10	40
6	20	30
9	30	20
12	40	10

3) Determination of specific gravity of experimental liquid:

Weight of empty sp.	Weight of sp. Gravity	Weight of sp. Gravity	Sp. Gravity of the
Gravity bottle (w ₁ g)	bottle+ water (w ₂ g)	bottle+ expt.	experimental solution
		Solution(w _x g)	$S_x = (w_x - w_1)/(w_2 - w_1)$
		w3=	$S_{3}=$
		$W_6 =$	$S_6 =$
		W9=	S9=
		$W_{12} =$	$S_{12} =$
		W ₁₅	$S_{15} =$

4) Determination of viscosity coefficient of experimental liquids:

Sl. No	Time of flow for fixed volume of liquid (sec.)						
	Water	Glycerol solution					
		3%	6%	9%	12%	15%	
(i) (ii)							
Average time of flow	$t_w =$	t ₃ =	t ₆ =	t ₉ =	t ₁₂ =	t ₁₅ =	

Calculation :

From the literature value of (η_w) at experimental temperature, absolute value of viscosity coefficient of experimental liquids can be determined,

$$\eta_x = \eta_w S_x \cdot t_L / t_w = [X = 3, 6, 9, 12, 15]$$

Graph Plotting

Coefficient of viscosity (poise)					
Concentration of Glycerol	3%	6%	9%	12%	15%
solution					

Plotting a graph of viscosity coefficient vs. concentration gives a calibration curve from which unknown concentration of glycerol solution can be obtained.

Experiment 2: Determination of partition coefficient for the distribution of I₂ between water and CCl₄.

Theory:

If a solute (soluble in both the solvents) is added to a pair of immiscible solvents then it will distribute itself between the two solvents in a constant temperature and pressure. This is Nernst distribution law. Let C and C be the molar concentration of the dissolved substance (solute) in solvent I and II respectively. Then according to the above law

$$\frac{C_1}{C_2} = K$$
, at constant temperature and pressure

Where K is a constant, called distribution or partition co-efficient of the dissolved substance between the solvent I and II.

The thermodynamic criterion for the equilibrium of the above system is that the chemical potential, μ_1 of the dissolved substance in solvent I must be equal to its chemical potential µ₂ in solvent II, i.e.,

$$\label{eq:main_state} \begin{split} \mu_1 = \mu_2 \\ or, \ \mu_1^*(T,P) + RTlna_1 = \mu_2^*(T,P) + RTlna_2 \end{split}$$

Where a_1 and a_2 are activities of the dissolved substance in solvent I and II respectively and μ_1^* & μ_2^* are constant temperature and pressure. Then it follows from the relation that $\frac{a_1}{a_2}$ = Constant, at constant temperature and pressure

If the two solutions behave ideally then the activities will be equal to the respective mole fraction (x), so that

 $\frac{a_1}{a_2} = \frac{x_1}{x_2} \approx \frac{c_1}{c_2}$, if the solutions are dilute. Then, $\frac{c_1}{c_2} = K$, at constant temperature and pressure provided the solutions are ideal and dilute and the dissolved substance exists in the same molecular form in the two solvents (and of course the two solvents must be immiscible, otherwise $\mu_1^* \& \mu_2^*$ will depend on the chemical nature of the solvent).

For example if iodine is added to a mixture of water and carbon tetrachloride and is shaken vigorously, at equilibrium, it will distribute itself between the two solvents in a constant ratio of concentration. Concentration of Iodine in each layer is determined by titrating against sodium thiosulphate solution using starch as indicator.

Apparatus required:

- 1) 250 ml stoppered glass bottles -2
- 2) 5 ml pipette -1
- 3) 25 ml pipette 1
- 4) Burette -1
- 5) 500 ml glass bottle -1
- 6) 250 ml volumetric flask 1
- 7) 250 ml conical flask -2

Chemicals required:

i)Saturated I₂ solution in CCl₄, ii) Pure CCl₄, iii) 10% KI, iv) 1% starch solution, v) Na₂S₂O₃.

Procedure:

- 1) Prepare 250 ml ~ (N/20) Na₂S₂O₃ solution.
- 2) Prepare 250 ml ~ (N/100) Na₂S₂O₃ solution by pipetting out 50 ml ~ (N/20) Na₂S₂O₃ solution in a 250 ml volumetric flask followed by dilution up to the mark with deionized water.
- 3) Prepare two sets in 250ml stoppered glass bottle as follows :

Set	Volume of saturated l solution in CCl4	2Volume of CCl4(ml)	Volume of water(ml)	Total volume(ml)
1	35	15	100	150
2	15	35	100	150

After mixing stopper the bottles properly and shake vigorously for an hour. Allow to settle for 10 minutes so that two layers become clearly separated.

- Pipette out 25 ml of the aqueous layer carefully from each set (so that no organic layer is taken out) and titrate with ~ (N/100) Na₂S₂O₃ solution using starch as indicator. Repeat one more time.
- 5) Pipette out 5 ml of the organic layer carefully from each set (so that no aqueous layer is taken out) in a 250 ml conical flask. Add 25ml water and 10 ml 10% KI solution to it. Shake well by swirling motion and titrate with ~ (N/20) Na₂S₂O₃ solution using starch as indicator. Repeat one more time.

Experimental result:

- 1) Room temperature:
- 2) Preparation of 250ml ~ (N/20) Na₂S₂O₃ solution

- 3) Prepare 250ml ~ (N/100) Na₂S₂O₃ solution from ~ (N/20) Na₂S₂O₃ solution:
- 4) Table 1:Titration of solvent layers

Set	Organic layer				Aqueous layer			
	No. of obs.	Vol. of organic layer(ml)	Vol. of (N/20) Na ₂ S ₂ O ₃ solution (ml)	~Mean vol.of ~ (N/20) Na ₂ S ₂ O ₃ solution (ml)	No. of obs.	Vol. of aqueous layer	Vol. of ~ (N/100) Na ₂ S ₂ O ₃ solution (ml)	Mean vol. of ~ (N/100) Na ₂ S ₂ O ₃ solution(ml)
1	i.	5			i.	25		
	ii.	5			ii.	25		
2	i.	5			i.	25		
	ii.	5			ii.	25		

Calculation:

If V₁ ml of 250ml ~ (N/20) Na₂S₂O₃ solution is required for 5 ml of organic layer and V₂ ml of ~ (N/100) Na₂S₂O₃ solution is required for 25 ml of aqueous layer then partition co-efficient K is given by

$$\mathbf{K} = \frac{C_1}{C_2} = \frac{\frac{V_1 \times (\frac{N}{20})}{5}}{\frac{V_2 \times (\frac{N}{100})}{5}} = 25(\frac{V_1}{V_2})$$

Set	Mean vol. of ~ (N/20) Na ₂ S ₂ O ₃ solution(ml) (V ₁)	Mean vol. of ~ (N/100) Na ₂ S ₂ O ₃ solution(ml)	$K=25(V_1/V_2)$	Mean K
I.		(V ₂)		
II.				-

Precautions and suggestions:

Care should be taken during pipetting out different layers so that no mixing occurs. Use bulb pipetizer instead of mouth sucking during pipetting out layers.

Experiment 3: Determination of Keq for $KI + I_2 = KI_3$, using partition coefficient between water and CCl_4

Theory:

In aqueous solution I2 and KI from complex as below

$$KI+I_2 = KI_2$$

For which equilibrium constant (Keq) is given by

$$K_{eq} = \frac{[KI_3]_{aq}}{[I_2]_{aq,free}[KI]_{aq,free}} \qquad \dots \dots \dots \dots (i) , \text{ for dilute solution}$$

This equilibrium constant can be determined by applying Nernst distribution law if the second immiscible organic solvent is so chosen that only one of the reactants or products can be distributed between two layers. In this case CCl_4 is chosen as second solvent in which I_2 is soluble. If distribution or partition coefficient (K_d)of I_2 in CCl_4 and water at a particular temperature is given by

$$\frac{[I_2]_{CCl_4}}{[I_2]_{water}} = K_a$$

Then at equilibrium, concentration of free iodine in aqueous layer,

$$[I_2]_{aq, \text{ free}} = \frac{[I_2]_{CCl_4}}{K_d}$$

[KI_3]_{aq} = Total I_2 in aqueous layer -[I_2]_{aq, \text{ free}}
= [I_2]_{aq, \text{total}} - \frac{[I_2]_{CCl_4}}{K_d}

Concentration of free KI in aqueous solution,

 $[KI]_{aq,free}$ = Total concentration of KI- $[KI_3]_{aq}$

$$= [KI]_{aq, total} - \{ [I_2]_{aq, total} - \frac{[I_2]_{CCl_4}}{K_d} \}$$

Therefore, $K_{eq} = \frac{[KI_3]_{aq}}{[I_2]_{aq, free}[KI]_{aq, free}} = \frac{[I_2]_{aq, total} - \frac{[I_2]_{CCl_4}}{K_d}}{\left(\frac{[I_2]_{CCl_4}}{K_d}\right)([KI]_{aq, total} - \{ [I2]_{aq, total} - \frac{[I_2]_{CCl_4}}{K_d} \})}$

Thus from the knowledge of $[I_2]_{aq, total}$, $[I_2]_{CCl4}$, $[KI]_{aq, total}$ and K_d , the value of K_{eq} can be determined.

Apparatus required:

- 1) 250 ml volumetric flask -2
- 2) Glass bottle -1
- 3) Burette -1

- 4) 25 ml pipette -1
- 5) 10 ml pipette -1
- 6) 250 ml stoppered glass bottle -2
- 7) 500 ml conical flask -1
- 8) 250 ml conical flask -2
- 9) Watch glass -1

Chemicals required: K₂Cr₂O₇, Na₂S₂O₃, KI, I₂ solution in CCl₄, 4(N) H₂SO₄, Starch. **Procedure:**

- Prepare 250ml of standard (N/20) K₂Cr₂O₇ solution and 250ml of standard (N/20) KI solution by accurate weighing.
- 2) Prepare 500ml ~ (N/20) Na₂S₂O₃ solution.
- 3) Standardize the (N/20) sodium thiosulphate solution idometrically using starch solution as indicator. Take 25ml of standard (N/20) K₂Cr₂O₇ solution in a 500 ml conical flask. Add about 10ml (one test tube) of 10% KI solution and 25ml of 4(N) H₂SO₄ solution. Cover the conical flask with watch glass and keep in dark for about 5 minutes. Add 150ml of distilled water and titrate the liberated iodine with sodium thiosulphate solution using starch as indicator.
- 4) Prepare two sets in dry, clean 250ml stoppered glass (leak proof) bottles:

Set	Volume of (N/20) KI solution(ml)	Volume of I ₂ solution in CCl ₄ (ml)	Volume of water(ml)
Ι	10	~40	90
II	20	~40	80

Stopper the bottles properly and shake the mixtures thoroughly for 45 minutes and allow to stand till complete separation of layers.

- 5) Pipette out 10 ml of the organic layer carefully from each set (so that no aqueous layer is taken out) in a 250ml conical flask. Add 25 ml water and 10ml 10% KI solution to it. Shaken well by swirling motion and titrate with standard (N/20) Na₂S₂O₃ solution using starch as indicator. Repeat one more time. Calculate concentration of I₂ in organic layer.
- 6) Pipette out 10ml of the aqueous layer carefully from each set (so that no organic layer is taken out), add ~ 40ml water and titrate with standard (N/20) Na₂S₂O₃ solution using starch as indicator. Repeat one more time. Calculate total concentration of I₂ in aqueous layer.
- 7) Calculate the value of K_{eq} using the supplied value of $K_{d.}$

Experimental result:

MIDNAPORE CITY COLLEGE

- 1) Room temperature.
- 2) Preparation of 500ml ~ (N/20) Na₂S₂O₃ solution.
- 3) Preparation of 250ml standard (N/20) K₂Cr₂O₇ solution.
- 4) Preparation of 250ml standard (N/20)KI solution.
- 5) Standardization of the prepared $Na_2S_2O_3$ solution:

Vol. of (N/20) K ₂ Cr ₂ O ₇ solution	Burette read	Burette reading		Average vol. of
(ml)	Initial(ml)	Final(ml)	$Na_2S_2O_3$ used	$Na_2S_2O_3$ (ml)
			(ml)	
25ml				
25ml				

Strength of Na₂S₂O_{3:}

6) Preparation of sets:

Set	Volume of (N/20) KI solution (ml)	Volume of I ₂ solution in CCl ₄ (ml)	Volume of water (ml)	[KI]aq,total

7) Titration of solvent layers:

Set	Organic Layer			Aqueous Layer				
	No. of Obs.	Vol. of organic layer (ml)	Vol. of standard (N/20) Na ₂ S ₂ O ₃ solution (ml)	Mean vol. of standard (N/20) Na ₂ S ₂ O ₃ solution [v _{org} (ml)]	No. of Obs.	Vol. of aqueous layer (ml)	Vol. of standard (N/20) Na ₂ S ₂ O ₃ solution (ml)	Mean vol. of standard (N/20) Na ₂ S ₂ O ₃ solution [v _{aq} (ml)]
Ι	i	10			i	10		
	ii	10]	ii	10		
II	i	10			i	10		
	ii	10]	ii	10]

Calculation:

For Set I :

Concentration of I₂ in organic layer, $[I_2]_{CCl_4} = \frac{Mean \ vol.of \ Na_2S_2O_3(V_{org}) \times its \ strength}{10}$ Total concentration of I₂ in aqueous layer, $[I_2]_{aq, \ total} = \frac{Mean \ vol.of \ Na_2S_2O_3(Vaq) \times its \ strength}{10}$ Similarly for set II calculate $[I_2]_{CCl_4}$ and $[I_2]_{aq, \ total}$. Calculate k_{eq} for set I and set II using the relation,

$$K_{eq} = \frac{[I_2]_{aq,total} - \frac{[I_2]_{CCl_4}}{K_d}}{\left(\frac{[I_2]_{CCl_4}}{K_d}\right) ([KI]_{aq,total} - \left\{[I_2]_{aq,total} \frac{[I_2]_{CCl_4}}{K_d}\right\}}$$

Experiment 4: Conductometric titration of an acid (strong, weak/ monobasic, dibasic) against base strong. Theory:

The conductance of an electrolyte solution depends on the mobility and concentration of each type of ions present in the solution. Any reaction which is associated either with replacement of one kind of ion by another of different mobility or with a charge I ionic concentration will be associated with a change in conductance of the reaction mixture. For such reaction gradual addition of one of the reactant to the other will result in a change of conductance in particular manner (depending on the nature of the reaction) up to equivalence point and in different manner after equivalence point. Then equivalence point corresponds to the break in titration curve obtained on plotting the conductance against the volume of the added substance (titrant).The conductance should be measured in a conductivity cell placed in a thermostat (to avoid change in conductance due to change in temperature). The titrant used should be many times stronger than the solution to be titrated to avoid change in conductance due to dilution .To avoid change in conductance due to change of electrolyte concentration resulting from electrolysis with direct current, measurement of conductance should be done using Alternating current based on Wheatstone Bridge principle. The electrodes used must be platinized (to avoid polarization). Different cases of Conductometric titration:

1) Strong acid vs. strong base:

Let us consider a strong acid (e.g.HCl) is titrated by a strong base (e.g. NaOH). Here the reaction involved is

$$H^++Cl^-+Na^++OH^- \rightarrow Na^++Cl^-+H_2O$$

In this reaction H⁺ ions are replaced by Na⁺ ions of lower mobility. Therefore conductance of the reaction mixture will gradually decrease linearly up to equivalence point and after equivalence point the conductance will increase linearly due to excess of added electrolyte (NaOH). The nature plot conductance vs volume of NaOH added is shown in the Figure-7.1, where intersection of two straight lines (say O) represents the equivalence point.



a) Titration of weak monobasic acid (CH₃COOH) by strong base (NaOH):

Here the reaction involved is

 $CH_{3}COOH \rightleftharpoons CH_{3}COO^{-} + H^{+} + Na^{+} + OH^{-} \rightarrow Na^{+} + CH_{3}COO^{-} + H_{2}O$

Initially H^+ ions (resulting from slight dissociation of CH_3COOH) are replaced by Na^+ ions of lower mobility. Due to this and also suppression of dissociation, conductance of reaction mixture initially slightly decreases. After appreciable addition of NaOH the conductance of reaction mixture will gradually increase up to equivalence point due to formation of highly ionized $CH_3COO^-Na^+$. After equivalence point the conductance will further increase at a different rate. The break in titration curve (Figure 7.2) indicates the equivalence point.



b) Titration of weak dibasic acid (Oxalic) by strong base (NaOH):

Here the reaction involved is



Oxalic acid is a dibasic acid and its $pk_2 > pk_1$. Therefore initially up to first equivalence point it will act as a relatively strong acid and conductance of the reaction mixture will gradually decrease linearly up to 1st equivalence point (O) (Figure 7.3) (because H⁺ ions are replaced by HC₂O₄. ions of lower mobility). Then conductance increases slowly up to 2nd equivalence point due to increase of $C_2O_4^{2-}$ ions with gradual addition of NaOH. After 2nd equivalence point (O₁) conductance will increase rapidly due to excess of added electrolyte(NaOH).



Apparatus required:

- 1) Conductivity meter
- 2) 100ml beaker
- 3) Burette (50ml), pipette (10ml) -1
- 4) 100ml volumetric flask -1
- 5) 250ml Glass bottle -3

Chemicals required:

Oxalic acid, NaOH, Acetic acid, HCl

Procedure:

- 1) Prepare 100ml \sim (N/2) NaOH solution.
- 2) Prepare 100ml (N/10) Oxalic acid solution by exact weighing in a volumetric flask.
- 3) Prepare 100ml ~(N/10) HCl and 100ml ~(N/10) acetic acid solution.
- 4) Take the NaOH solution in a burette. Determine the volume of 50 drops and calculate the volume of one drop.
- 5) Standardize the \sim (N/2) NaOH solution in the following way. Pipette out 10ml of the prepared oxalic acid solution in a 100ml beaker, add sufficient distilled water

to cover the electrodes of the conductivity cell. Stir the solution gently and take the conductance reading. Add (N/2) NaOH solution from a burette, two drops at a time initially and subsequently one drop at a time. Stir gently and measure conductance after each addition. After end point take at least 6-8 readings.

- 6) Draw the graph of conductance vs number of drops of NaOH added . From the second equivalence point calculate strength of NaOH solution.
- 7) Standardize the prepared HCl and acetic acid solution following the procedure of step 5.
- 8) From the graph of conductance vs number of drops of standard NaOH added, calculate strength of prepared HCl and acetic acid solution.

Experimental result:

- 1) Room temperature.
- 2) Preparation of 100ml 0.1(N) Oxalic acid.
- 3) Preparation of 100ml approximately (N/2) NaOH solution.
- 4) 50 drops of NaOH solution=..... ml NaOH.
- 5) Conductometric titration: Oxalic acid vs NaOH:

Volume of (N/10) oxalic acid	No. of drops of NaOH added	Total no, of drops of NaOH added	Conductance (mho)
10ml			

Total no. of drops of NaOH corresponding to second equivalence point \equivml NaOH(V₁).

Strength of NaOH solution = $\frac{10 \times 0.1(N)}{\text{Volume of NaOH}(V1)}$ 6) Conductometric titration : HCl vs NaOH

Calculate strength of HCl in a similar way as step 5.

Volun	ne of HCl	No. NaC	of drops of H added	Total no. of drops of NaOH (Conductance(mho)
		1 140	11 uuuuu	aaaee	•	
	Volume of ace	tic	No. of drops of		Total no. of drops of	Conductance(mho)
	acid		NaOH added		NaOH added	
	10ml					

7) Conductometric titration: Acetic acid vs NaOH

Calculate strength of Acetic acid in a similar way as step 5

Precautions and Suggestions:

- 1) Add NaOH solution from a burette, two drops at a time initially and subsequently one drop at a time.
- 2) Stir gently before taking each conductance reading.

Experiment 5: Study of saponification reaction conductometrically Theory:

Saponification of ester refers to the aqueous alkaline hydrolysis of any type of ester. As an example saponification of methyl acetate can be represented as

 $CH_3COOCH_3 + OH^- \rightarrow CH_3COO^- + CH_3OH$

The rate of the reaction, $R = \frac{d[ester]}{dt} = k[ester][OH^-]$, i.e.

the overall kinetic order of the reaction is two.

If 'a' be the initial concentration of ester and 'x' be the concentration of ester reacted at time t then the concentration of unreacted ester is (a-x). For equal initial concentration of both the reactants we have,

 $\frac{d(a-x)}{dt} = k(a-x)^2$, Where k is the second order rate constant with unit of mol⁻¹L.s⁻¹ Or, $\frac{dx}{dt} = k(a-x)^2$

Integrating this equation by imposing boundary condition, at t=0, x=0, we have,

$$\frac{x}{(a-x)} = k$$
 at

The course of the reaction can be monitored by measuring conductance of the reaction mixture with time. In this case the conductance of the reaction mixture decreases with time as the OH⁻ ions are replaced by CH₃COO⁻ ions of lower mobility. Therefore from the measured conductance value the relative amounts of the ions present and hence the extent of reaction can be predicted. If C₀, C_t, and C ∞ are the values of conductance at time t=0, t=t and t= ∞ respectively, then

$$a \propto (C \propto - C_0)$$

(a-x) \le (C \le - C_t)
x \le (C_t - C_0)
get,
$$\frac{C_0 - C_t}{C_t - C} = k \text{ at}$$

Substituting in equation (i) we get, $\frac{c_0 - c_t}{c_t - c_\infty} =$

Plot of $\frac{C_0 - C_t}{C_t - C_{\infty}}$ vs t should yield a straight line passing through the origin with a slope 'ka'. From the known value of 'a', the rate constant for the saponification reaction, k can be calculated.

Apparatus required:

1) 100ml vol. flask -5

- 2) 250ml conical flask -1
- 3) 100ml beaker -1
- 4) 500ml glass bottle -1
- 5) Burette -1
- 6) Pipette 25ml -2
- 7) Pipette 10ml -1
- 8) Pipette 1ml -1

Chemicals required:

Oxalic acid, NaOH, Acetic acid, Phenolphthalein, Methyl acetate

Procedure:

- 1) Prepare 100ml of \sim (N/10) oxalic acid solution by exact weighing.
- Prepare 250ml of ~(N/10) NaOH solution. Standardize this NaOH solution against 10ml standard oxalic acid solution.(Use Phenolphthalein).
- 3) Dilute the standardized NaOH solution to prepare exact (M/60)NaOH solution of volume 100ml using deionized water.
- Prepare 100ml of a ~(N/10) acetic acid (HAc) solution. Standardized 10ml of it against the standardized ~(N/10) NaOH solution (Use Phenolphthalein). From the standardized solution prepare 100ml exact (M/60) HAc solution by proper dilution.
- 5) Take a 100ml volumetric flask with about 50ml of deionized water. Transfer into this exactly 1.0ml of pure methyl acetate, using a 1ml graduated pipette, and make up the volume with deionized water. Estimate the concentration of the solution in molarity using the following data:

Molecular weight of methyl acetate = 74.08Density of the supplied methyl acetate at room temperature (t⁰C)

=0.932- (t-20) × 1.25×10^{-4} g/ml

Prepare 100ml exact (M/60)methyl acetate solution from this standard methyl acetate solution by proper dilution with deionized water.

- Prepare 50ml exact (M/120) NaOH solution from exact (M/60) NaOH solution by proper dilution (use a 25ml pipette) with deionized water and note its conductance as C₀.
- 7) Mix 25ml each of (M/60) NaOH and (M/60) HAc solutions and note its conductance as $C\infty$.
- 8) Take 25ml (M/60) methyl acetate and add 25ml (M/60) NaOH solution to it. Note the time of half – discharge of the pipette and homogenize the mixture. Measure the conductance (Ct) of the solution at intervals of about 1 minute. Take at least 15 reading.

9) Plot $(C_0-C_t)/(C_t-C_\infty)$ versus time to obtain the rate constant of the reaction at room temperature.

Experimental result:

- 1) Room temperature.
- 2) Preparation of 100ml 0.1 (N) oxalic acid.
- 3) Preparation of 250ml approximately 0.1(N) NaOH solution.
- 4) Standardisation of NaOH solution against standard oxalic acid solution.

Volume of oxalic acid (ml)	Burette reading		Volume of NaOH used	Average volume of NaOH (ml)
	Initial(ml)	Final(ml)		
10				
10				

Strength of NaOH solution.

- 5) Preparation of 100ml exact (M/60) NaOH solution.
- 6) Preparation of 100ml ~ $\left(\frac{N}{10}\right)$ acetic acid solution.
- 7) Standardisation of HAc solution against standard NaOH solution.

Volume of HAc (ml)	Burett	e reading	Volume of NaOH	Average volume of NaOH	
	Initial(ml)	Final(ml)	used (ml)	(ml)	
10					
10					

Strength of HAc solution.

- 8) Preparation of 100ml exact (M/60) HAc solution.
- 9) Preparation of 100ml exact (M/60) methyl acetate solution.
- 10) Measurement of C_0 and $C\infty$

C₀=conductance of (M/120) NaOH solution.

 $C\infty$ = conductance of (M/120) sodium acetate solution.

11)Recording of conductance data for saponification reaction.

Conclusion:

Time (t) in minute	Conductance	$(C_t)C_0-C_t(mho)$	C_t - $C\infty$ (mho)	$(C_0-C_t)/(C_t-C\infty)$
	(mho)			

Plot of $\frac{c_0 - c_t}{c_t - c_{\infty}}$ vs t gives a straight line passing through the origin. From the graph

slope = ka.

Or, k=Slope/a, where 'a' is initial concentration of reactants = (M/120).

Experiment 6: Verification of Ostwald's dilution law and determination of Ka of weak acid.

Theory:

Electrolytes obey Ohm's law i.e.,

Current (I) \propto Voltage (V) applied

The proportionality constant is called conductance (C) [inverse of resistance (R)] of the conductor concerned. The conductance of an electrolyte solution is measured with the help of a conductivity cell. If two parallel electrodes, each of area 'A' are placed 'I' distance apart then the measured conductance (C) of an electrolyte solution will be

 $C = \frac{1}{R} = \frac{1}{\rho} \frac{A}{l}$, ρ is the resistivity

 $= k \frac{A}{i}$, k is the specific conductance or conductivity of the electrolyte solution.

Or, $k = C\frac{l}{A}$, Where $\frac{l}{A}$ is called cell constant of the conductivity cell

i.e., Specific conductance = Conductance \times Cell constant

The specific conductance (k) of an electrolyte solution is its conductance when unit volume of the solution is placed between two parallel electrodes each of unit area and set unit distance apart. In C.G.S. its unit is $ohm^{-1} cm^{-1}$.

Cell constant of a conductivity cell is determined by measuring the conductance of an electrolyte solution of known specific conductance at the experimental temperature within the same conductivity cell.

Equivalent conductance (λ) of an electrolyte solution is its conductance associated with a definite volume of the solution containing 1 gm equivalent of the electrolyte placed between two large electrodes set unit distance apart. Therefore,

Equivalent conductance (λ) = Specific conductance (k) × volume of the solution in cm³ containing 1 gm equivalent of the electrolyte

Unit of ' λ ' in CGS in ohm⁻¹cm²g equi⁻¹.

Equivalent conductance of a weak electrolyte will increase with dilution and this is mainly due to increase n degree of dissociation, the variation of ionic mobility with dilution being insignificant due to low concentration of the ions.

Let α be the degree of ionization of a weak acid HA at concentration C equivalent/lit.

Then at equilibrium

$$HA \rightleftharpoons H^+ + A^-$$

C(1-\alpha) C\alpha C\alpha

The dissociation constant or ionization constant (Ka) can be represented as

$$K_a = \frac{C^2 \alpha^2}{C(1-\alpha)} = \frac{C \alpha^2}{1-\alpha}$$

If equivalent conductance of the weak electrolyte at concentration/ lit and at infinite dilution are λ and λ_0 respectively then,

$$\alpha = \frac{\lambda}{\lambda_0}$$

so, $K_a = \frac{C(\lambda/\lambda_0)^2}{1 - \frac{\lambda}{\lambda_0}}$
or, $1 - \frac{\lambda}{\lambda_0} = \frac{C\lambda^2}{K_a \lambda_0^2}$ or, $\frac{1}{\lambda} = \frac{1}{\lambda_0} + \frac{\lambda C}{\lambda_0^2 K_a}$

This is known as Ostwald dilution law. The plot of $\frac{1}{\lambda}$ vs. λC will be a straight line with a positive intercept $\frac{1}{\lambda_0}$ from which λ_0 can be calculated. From the slope $\frac{1}{\lambda_0^2 K_a}$, K_a can be calculated using the temperature corrected literature value of λ_0 .

Apparatus required:

- 1) 100ml vol. flask -3
- 2) 250ml conical flask -1
- 3) 100ml beaker -2
- 4) 500ml glass bottle -2
- 5) Burette -1
- 6) Pipette 25 ml -2
- 7) Pipette 10ml -1
- 8) 250ml vol. flask -2

Chemicals required:

Oxalic acid, Acetic acid, NaOH, KCl, Phenolphthalein

Procedure:

- 1) Prepare 100ml \sim (N/10) oxalic acid solution by accurate weighing.
- 2) Prepare approximately 250ml of \sim (N/10) NaOH solution and standardize the NaOH solution against the prepared standard oxalic acid solution.(use phenolphthalein).
- 3) Prepare 250ml of \sim (N/10) acetic acid solution in conductivity water and standardize the solution against the standard NaOH solution (use phenolphthalein).
- 4) Prepare 250ml ~ 0.1 (N) [slightly higher than 0.1 (N) KCl solution by accurate weighing and prepare 100ml of an exact 0.1 (N) KCl solution by proper dilution. Prepare 100ml of an exact 0.01 (N) KCl from the exact 0.1 (N) KCl solution. Determine the cell constant of the conductivity cell using the exact 0.1 (N) KCl and the prepared exact 0.01 (N) KCl solution .With the help of the literature value of specific conductances of these solutions at room temperature,

calculate the mean value of cell constant and use it subsequently. Measure the conductance of conductivity water also.

- 5) Prepare 250ml, exact (N/50) weak acid solution from the standardized solution using conductivity water. With the help of 25ml pipette, take 50ml of this solution in the clean and dry conductivity cell and measure its conductance.
- 6) Use the same 25ml pipette to take out 25ml of the (N/50) weak acid solution from the conductivity cell. Pipette out 25ml of the conductivity water into the conductivity cell to make the solution exactly (N/100) in situ. Mix the solution well by careful swirling (so that no solution comes out). Measure the conductance and note it as that of exact (N/100) weak acid solution. Separate pipettes for weak acid solutions and conductivity water may be used.
- 7) Follow the procedure of step (6) to prepare in situ exact (N/200), (N/400), (N/800), and (N/1600) weak acid solutions in steps and note their conductances.
- 8) Calculate the equivalent conductivities of the diluted solutions of the weak acid using the mean value of cell constant. Apply corrections for specific conductance of conductivity water.
- 9) From the plot of $1/\lambda$ versus (λC) calculate the λ_0 from the intercept and from the temperature corrected literature value of ion conductances. Calculate K_a and pk_a of the weak acid from the slope using the temperature corrected literature value of λ_0 .

Experimental result:

- 1) Room temperature.
- 2) Preparation of 100ml 0.1(N) oxalic acid.
- 3) Preparation of 250ml approximately 0.1(N)NaOH solution.
- 4) Standardisation of NaOH solution against standard oxalic acid solution:

Burette reading		Volume of NaOH	Average volume	
Initial(ml)	Final(ml)	used	of NaOH (ml)	
	Burette Initial(ml)	Burette reading Initial(ml) Final(ml)	Burette reading Volume of NaOH Initial(ml) Final(ml)	

Strength of NaOH solution:

- 5) Preparation of 250ml \sim (N/10) acetic acid solution:
- 6) Standardisation of acetic acid solution against standard NaOH solution :

Strength of acetic acid solution:

7) Preparation of 250ml exact (N/50) acetic acid solution.

- 8) Preparation of 250ml \sim 0.1 (N) [slightly higher than 0.1 (N) KCl solution.
- 9) Preparation of 100ml exact 0.1 (N) KCl solution and preparation of exact 0.01 (N) KCl .

Volume of acetic acid	Burette	e reading	Volume of NaOH	Average volume of NaOH (ml)
solution(ml)	Initial(ml)	Final(ml)	used	
10				
10				

10) Determination of cell constant:

Concentration of KCl solution	Conductance (mho)	Specific conductance (mho.cm ⁻¹)	Cell constant (cm ⁻¹)	Mean cell constant (cm ⁻¹)
0.01(N)				
0.1(N)				

11) Determination of equivalent conductance of different acetic acid solutions:

Concentration of	Observed	Corrected specific	Equivalent	1/ λ (mho ⁻	$\Lambda c(mhocm^{-1})$
acetic acid	conductance	conductance	conductance	1 cm $^{-2}$ eq)	
solution	(mho)	$(mho.cm^{-1})$	(λ) (mhocm ² eq ⁻		
		{Observed	1)		
		conductance-			
		Conductance of			
		conductivity			
		water) × mean			
		cell constant}			

The plot of $\frac{1}{\lambda}$ vs λ C will be a straight line with a positive intercept $\frac{1}{\lambda_0}$ from which λ_0 can be calculated.

From the slope $\frac{1}{\lambda_0^2 K_a}$, K_a can be calculated using the temperature corrected literature value of λ_0 . Temperature corrected literature value of $\lambda_0 = \lambda_t^0$ (H⁺) + λ_t^0 (OAc⁻) $\lambda_t^0 (\mathrm{H}^+) = \lambda_{25}^o [1 + 1.42 \times 10^{-2} (\mathrm{t-25})]$ $\lambda_t^0 (\mathrm{OAc}^-) = \lambda_{25}^o [1 + 0.02 (\mathrm{t-25})]$

C6P: Inorganic Chemistry Lab

Experiment 1: Estimation of available chlorine in bleaching powder.

Aim

To determine the available chlorine in the given sample of bleaching powder by the iodometric method.

Theory:

Bleaching powder is used as a bleaching agent and also as a disinfectant. The main constituent of bleaching powder is calcium hypochlorite $[Ca(OCl)_2]$ which supplies chlorine $[Cl_2]$ with dilute acids. So, the available chlorine is defined as the percentage of chlorine made available by bleaching powder when treated with dilute acids. The available chlorine present in bleaching powder sample is determined iodometrically by treating its solution with an excess of potassium iodide solution in the

acidic medium. The liberated iodine (I_2) is treated with sodium thiosulphate $(Na_2S_2O_3)$ solution using freshly prepared starch solution as indicator to be added near the end point.

Apparatus:

Digital Balance, Burette, Conical flask, Measuring flask, Funnel, Glass rod, Beakers Reagents:

Bleaching powder, Standard sodium thiosulphate solution

10% Potassium iodide (KI), dilute acetic acid, freshly prepared starch solution as indicator Bleaching powder is commonly used as a disinfectant. The chlorine present in the bleaching powder gets reduced with time. So, to find the exact quantity of bleaching powder required, the amount of available chlorine in the sample must be found out.

Chlorine will liberate free iodine from potassium iodide solution when its pH is 8 or less. The iodine liberated, which is equivalent to the amount of active chlorine, is titrated with standard sodium thiosulphate solution using starch as indicator.

Apparatus:

Mortar and pestle

- Volumetric flask
- > Burette

Reagents: Concentrated glacial acetic acid Potassium iodide Starch indicator

- Erlenmeyer flask.
- > Pipette

Procedure

Dissolve 1g bleaching powder in 1 litre of distilled water in a volumetric flask, and stopper the container.

(This can be done by first making a paste of the bleaching powder with mortar and pestle.)

Place 5 mL acetic acid in an Erlenmeyer flask and add about 1g potassium iodide crystals. Pour 25 mL of bleaching powder solution prepared above and mix with a stirring rod.

Titrate with 0.025 N sodium thiosulphate solution until a pale yellow colour is obtained. (Deep yellow changes to pale yellow.)

Add 1mL of starch solution and titrate until the blue colour disappears.

Note down the volume of sodium thiosulphate solution added (V_1) .

Take a volume of distilled water corresponding to the sample used.

Add 5 mL acetic acid, 1g potassium iodide and 1 mL starch solution.

If blue colour occurs, titrate with 0.025 N sodium thiosulphate solution until the blue colour disappears.

Record the volume of sodium thiosulphate solution added (A₁).

If no blue colour occurs, titrate with 0.025 N iodine solution until a blue colour appears. Note down the volume of iodine (A_2) .

Then, titrate with 0.025 N sodium thiosulphate solution till the blue colour disappears. Record the volume of sodium thiosulphate solution added (A₃). Note down the difference between the volume of iodine solution and sodium thiosulphate as A_4 ($A_4=A_2-A_3$).

Note: Blank titration is necessary to take care of the oxidising or reducing reagents' impurities.

Observation

Bleaching powder solution x Standard sodium thiosulphate solution (0.025 N)

Introduction

In acid solution practically all oxidizing agents will oxidize iodide ion to iodine quantitatively. The iodine formed in the reaction can then be titrated by means of a standard sodium thiosulfate solution. This type of indirect titration is given the general name of **iodometry**.

Iodometric methods of analysis have a wide applicability for the following reasons

1. Potassium iodide, KI, is readily available in high purity.

2. A good indicator, starch, is available to signal the equivalence point in the reaction between iodine and thiosulfate. Starch turns blue-black in the presence of iodine. Therefore, when the blue-black color disappears, the iodine has been completely reduced to the iodide ion.

3. Iodometric reactions are rapid and quantitative.

4. A precise and stable reducing agent, sodium thiosulfate ($Na_2S_2O_3$), is available to react with the iodine.

The amount of iodine liberated in the reaction between iodide ion and an oxidizing agent is a measure of the quantity of oxidizing agent originally present in the solution. The amount of standard sodium thiosulfate solution required to titrate the liberated iodine is then equivalent to the

Volume of bleaching powder	Durette	reading	Volume of titrant (mL)
solution (Inc)	Initial	Final	
			0
			5
			ц
			μ
			Initial Final

Distilled water \times Standard sodium thiosulphate solution (0.025 N)

Trial no.	Volume of distilled water (mL)	Burett	e reading	Volume of titrant (mL)
		Initial	Final	
	8			
	0			

Distilled water x Standard iodine solution (0.025N)

amount of oxidizing agent. Iodometric methods can be used for the quantitative determination of strong oxidizing agents such as potassium dichromate, permanganate, hydrogen peroxide, cupric ion and oxygen.

As has been mentioned above, the endpoint in a titration of iodine with thiosulfate is signaled by the color change of the starch indicator. When starch is heated in water, various decomposition products are formed, among which is beta-amylose which forms a deep blue-black complex with iodine. The sensitivity of the indicator is increased by the presence of iodide ion in solution.

Trial no.	Volume of distilled water (mL)	Burett	e reading	Volume of titrant (mL)
1		Initial	Final	
		2.		
		<i>z.</i>		

Calculation

 $(V - A_1)$ or $(V + A_4) \ge N \ge 35.46$

mg of $Cl_2/mL(B) =$

mL of bleaching powder solution taken

1000 mL of bleaching powder solution contains 1000 x B mg of Cl_2 i.e., 1000 mg bleaching powder contains 1000 B mg of Cl_2

Therefore, 100 mg of bleaching powder contains = $\frac{1000 \text{ x B}}{10}$

% of chlorine available =

Results

Available chlorine in the given bleaching powder is. ...%

Experiment 2: Iodometric Determination of Cu in Brass

However, if the starch indicator solution is added in the presence of a high concentration of iodine, the disappearance of the blue-black color is very gradual. For use in indirect methods, the indicator is therefore added at a point when virtually all of the iodine has been reduced to iodide ion, causing the disappearance of the color to be more rapid and sudden. The starch indicator solution must be freshly prepared since it will decompose and its sensitivity is decreased. However, a properly prepared solution will keep for a period of a few weeks. A preservative such as a small amount of mercuric ions may be added to inhibit the decomposition.

Solutions of sodium thiosulfate are made up to an approximate concentration by dissolving the sodium salt in water that has previously been boiled. Boiling the water is necessary to destroy

micro-organisms which metabolize the thiosulfate ion. A small amount of Na_2CO_3 is added to the solution in order to bring the pH to about 9. The solution is standardized by taking a known amount of oxidizing agent, treating it with excess iodide ion and then titrating the liberated iodine with the solution to be standardized. Oxidizing agents such as potassium dichromate, bromate, iodate or cupric ion can be employed for this procedure. You will be using potassium iodate, KIO₃, as your primary standard. The reaction between IO_3^- and I^- is given as

 $6H^++IO_3^-+5I^---->3I_2+3H_2O$

Reactions Involved in Iodometric Processes

Iodometric methods depend on the following equilibrium:

 $I_2 + I^{-} <==> I_3^{-}$

Since the solubility of I_2 in water is quite low, the formation of the tri-iodide ion, I_3 , allows us to obtain useful concentrations of I_2 in aqueous solutions. The equilibrium constant for this reaction is approximately 700. For this reason iodometric methods are carried out in the presence of excess iodide ion.

The reaction between iodine and the thiosulfate ion is:

 $I_2 + 2S_2O_3^{2-} <==> 2I^- + S_4O_6^{2-}$

This reaction proceeds quantitatively in neutral or slightly acidic solutions. In strongly alkaline or acidic solutions the oxidation of the thiosulfate does not proceed by a single reaction. In the former, the thiosulfate ion is oxidized to sulfate as well as to the tetrathionate. In the latter, the thiosulfuric acid formed undergoes an internal oxidation-reduction reaction to sulfurous acid and sulfur. Both of these reactions lead to errors since the stoichiometry of the reactions differs from that shown above for the thiosulfate as a reducing agent. The control of pH is clearly important. In many cases the liberated iodine is titrated in the mildly acidic solution employed for the reaction of a strong oxidizing agent and iodide ion. In these cases the titration of the liberated iodine must be completed quickly in order to eliminate undue exposure to the atmosphere since an acid medium constitutes an optimum condition for atmospheric oxidation of the excess iodide ion.

The basic reaction in the determination of copper using the iodometric method is represented by the equation:

 $2Cu^{2+} + 4I^{-} \ll 2CuI(s) + I_{2}$

This is a rapid, quantitative reaction in slightly acidic solutions, if there is a large excess of iodide ion present and if the copper is in the form of a simple ion rather than a complex one. The iodine that is liberated can be titrated in the usual manner with standard thiosulfate solution. The reaction involving cupric ion and iodide takes place quantitatively since the cuprous ion formed as result of the reduction is removed from the solution as a precipitate of cuprous iodide.

Iron interferes since iron(III) ions will oxidize iodide. Since the iron will be found in the +3 oxidation state as a result of the dissolution of the brass sample, a means of preventing this interference is necessary. This can be accomplished by converting the iron(III) to a soluble iron(III) phosphate complex using phosphoric acid. At a pH of 3.0-4.0 the iron phosphate complex is not reduced by iodide ion. If arsenic and antimony are present, they will provide no interference at this pH if they are in their higher oxidation states. Brass formulations also may contain up to 39% Zn,

2.5% Sn and 8.5% Pb. When dissolved in concentrated nitric acid, the zinc and the lead become Pb^{2+} and Zn^{2+} . These do not interfere with the analysis of copper because they are not reduced to the Pb^+ and Zn^+ states by the action of iodide ion under the conditions of this experiment. The tin is oxidized to Sn^{4+} by the concentrated nitric acid and after dilution and adjustment of pH this form becomes SnO_2 which is insoluble and may be observed as an inert white precipitate at the bottom of your flask. Under these conditions the tin does not interfere with the analysis.

Sources of Error

The following are the most important sources of error in the iodometric method:

1. Loss of iodine by evaporation from the solution. This can be minimized by having a large excess of iodide in order to keep the iodine tied up as tri-iodide ion. It should also be apparent that the titrations involving iodine must be made in cold solutions in order to minimize loss through evaporation.

2. Atmospheric oxidation of iodide ion in acidic solution. In acid solution, prompt titration of the liberated iodine is necessary in order to prevent oxidation.

3. Starch solutions that are no longer fresh or improperly prepared. The indicator will then not behave properly at the endpoint and a quantitative determination is not possible.

EXPERIMENTAL

Preparation of a 0.10 M Standard Na₂S₂O₃ Solution

With a graduated cylinder measure out 1 liter of distilled water. Place it in your 1 liter beaker and boil the water for at least 5 minutes. Weigh out 25 g of $Na_2S_2O_3 \cdot 5H_2O$ and 0.1 g of Na_2CO_3 . Dissolve the thiosulfate in the hot water and then cool this solution with the aid of an ice bath to room temperature. Then add the carbonate and stir until it is completely dissolved. Transfer the solution to your plastic 1 liter bottle. When not in use store this bottle in the darkness of your equipment cabinet as the decomposition of thiosulfate is catalyzed by light.

Blank Determination

Potassium iodide may contain appreciable amounts of iodate ion which in acid solution will react with iodide and yield iodine. The liberated iodine would react with thiosulfate and thereby cause the apparent molarity of the thiosulfate to be too low. The following procedure allows for the determination of a blank correction which will properly correct for any iodate that might be present. Prepare a solution of exactly 2.00 g of KI dissolved in 50 mL of distilled water and then acidify the solution with 5 mL of 3 M sulfuric acid and then immediately add 5 mL of starch indicator. If a blue-black color appears right after mixing, use the thiosulfate solution in the buret to determine the volume of solution required to cause the color- to disappear. This volume must be subtracted from the standardization and analyses volumes. If the potassium iodide is completely iodate-free no color will of course develop and no blank correction is necessary.

Standardization of the Na₂S₂O₃ Solution

Dry approximately 2 g of potassium iodate, KIO₃, at a temperature of 110 °C for one hour. Weigh to a precision of ± 0.0001 g three samples of the potassium iodate having weights near 0.12 g directly into three 250 mL Erlenmeyer flasks. Dissolve the iodate in 75 mL of distilled water. Cover the flasks with parafilm and store them. Rinse and fill your buret with the solution. Add 2.00 g of KI to each of the potassium iodate solutions. If a blank correction is required add exactly 2.00 g of KI to each. If no blank determination is required, the exact amount of KI is not crucial but should be close to 2 g. Then add 10 mL of 1 M HCl to one of the solutions. It will turn a darkbrown color. Immediately titrate it with the thiosulfate solution. When the color of the solution becomes very pale yellow add 5 mL of starch indicator. Continue the titration until the blue color of the starch complex just disappears. Follow the same procedure with each of the other two solutions, first adding the HCl then titrating. Correct your titration data for buret error and if necessary apply the blank correction. Calculate the molarity of the Na₂S₂O₃ solution. Results should agree to within 0.2% of the average. If you do not achieve that kind of precision, titrate additional samples.

Dissolution of the Brass Sample

The following procedures in this section make use of the hot plates in the fume hoods. The solutions of dissolved brass generally have a low volume and high acid and salt concentrations. "Bumping" or little explosions of steam in the superheated liquid can occur. You don't want your hand to be close to the mouth of the flask should the solution suddenly "bump" because drops of acid (not to mention part of your sample) will fly out of the flask and possibly onto your hand. For that reason you must use your tongs to place the flasks on the hot plate and to remove them. Don't use strips of paper towel or the rubber Hot Hands because your real hand will end up being too close to the mouth of the flask.

The brass sample which you will receive does not have to be dried before use. Accurately weigh out three brass samples, of about 0.3 g each, directly into separate 250 mL Erlenmeyer flasks. In the fume hood add 5 mL of 6 M HNO₃. Warm the solution on a hot plate in the fume hood until dissolution is complete. Add 10 mL of concentrated (not 3 M) H₂SO₄ and continue heating until white SO₃ fumes appear. It is important that you do not mistake ordinary water vapor for SO₃ fumes. It is also important at this point that the flask *not* be removed from the hood. SO₃ fumes are dangerous and ought not to be inhaled. Only when the slightly denser white fumes of SO₃ are observed can you be sure that all HNO₃ has been removed. NO₃⁻ will oxidize I⁻ and hence will seriously interfere with the procedure. Cool the flask in air for one or two minutes and then in an ice bath, then carefully add 20 mL of distilled H₂O. Boil for one or two minutes then again cool in an ice bath. Continue to keep the flask in the ice bath and using your medicine dropper add

concentrated NH₃(aq) dropwise, and with adequate mixing, until the light-blue color of the solution is completely changed to the dark-blue color of the copper tetraammine complex. As many as 400 drops (20 mL) may be required. The solution must be kept cool in an ice bath since the reaction between the concentrated H₂SO₄ and concentrated NH₃ is highly exothermic. Now add 3 M H₂SO₄ dropwise until the dark-blue color just disappears. You don't have to produce a complete disappearance of the dark blue color but you need to approach that point. The subsequent addition of phosphoric acid will lower the pH appropriately to around 3.5. If you add too much 3M H₂SO₄ the pH may turn out to be sufficiently low to cause unwanted side reactions to occur when you reduce the Cu²⁺ with iodide. If you are uncertain about the disappearance of the dark blue color you may put 50 mL of 0.06 M Cu²⁺ in a clean 250 mL flask and add 12 M ammonia dropwise until you have that unmistakable dark blue color. Then add 3M H₂SO₄ dropwise until the blue color *almost* disappears. Then add 2 mL concentrated phosphoric acid and you ought to see the dark color completely disappear. You may copy that procedure to achieve an appropriate pH of around 3.5 for subsequent steps in the analysis. Now, back to your real sample: Once you are confident that you haven't added too much 3M H₂SO₄, but that you have caused the dark color of the copper tetraammine complex almost to disappear, add 2.0 mL of concentrated phosphoric acid, H₃PO₄, to each sample. Verify to yourself that they exhibit the light copper color rather than the dark color and cover the flasks with parafilm and set them aside until you are ready to proceed with the titration.

Titration of the Dissolved Brass Sample

If you have let the dissolved samples stand overnight, be sure to warm the sample on a hot plate (this can be done at your desk) in order to dissolve all larger crystals of copper sulfate that might have formed. Be sure to cool the samples to room temperature, or below, with the aid of an ice bath. The solutions will still contain a fine, white precipitate at this point; however, this will not interfere with the rest of the procedure. From this point on work with only one sample at a time. Add 4.0 g of KI to one of your samples and titrate immediately with the standard thiosulfate solution. The sample contains white CuI precipitate and the color of I_3^- must be observed against that precipitate. The slurry will at first appear brown or dark yellow-brown. Continue adding thiosulfate until the slurry is a light mustard color. At this point add 5 mL of starch indicator and titrate until the mixture in the flask takes on a milky pink or lavender hue. Now add 2 g of KSCN and mix well; the solution will darken somewhat. After the addition of thiocyanate, continue to add more thiosulfate dropwise. You should observe a sudden change to a white or cream color. That is the endpoint of the titration. After you have titrated all three samples calculate the percentage of Cu in each of the brass samples, the average percentage and the average deviation. The description above applies for brass samples with low concentrations of zinc (<10%). Some of you may have brass samples with higher concentrations of zinc. Such samples will become quite dark after the addition of KI and will lighten only slightly as thiosulfate is added. The "mustard color" will be darker than samples having low percentages of copper. When the starch is added the sample will become dark blue-black again and as you approach the end point with the thiosulfate the slurry will turn a violet color rather than milky pink or lavender hue. With the addition of the

KSCN the solution will darken somewhat as in the case of the other samples, but the final end point will be a bit darker than the white or cream color described above. If you think that you have a sample with high zinc content, observe your progress carefully and take notes which will allow you to achieve repeatability.

Explanation: The reduction of Cu^{2+} to Cu^+ occurs as the result of the oxidation of I⁻ to I₂. The I₂ combines with iodide ion to produce the dark brown triodide ion, I₃⁻. The excess iodide ion also causes the reduced copper to precipitate as white cuprous iodide, CuI. I₂ and I₃⁻ in solution tend to adsorb on the surface of CuI thus becoming unavailable for rapid reduction by the thiosulfate. As a result, iodometric titrations involving reduced copper tend to yield lower results unless the adsorbed I₂ can be liberated by adding thiocyanate ion, SCN⁻, which competes with the adsorbed iodine molecules on the surface of solid particles of CuI. After the addition of thiocyanate, continue to add more thiosulfate dropwise. You should observe a sudden change to a white or cream color. That is the endpoint of the brass samples, the average percentage and the average deviation. REPORT

Your report must include the following information in the two sections below.

- 1. Unknown number
- 2. The three weights of KIO3 used for the standardization of thiosulfate
- 3. Volume in mL of thiosulfate for each of your standardization titrations
- 4. Average molarity of the thiosulfate solution
- 5. Weight of brass used for each sample
- 6. Volume of thiosulfate solution used for each sample
- 7. Cu percentage for each sample
- 8. The average percentage of Cu in your brass sample
- 9. The average deviation from the mean of the percent Cu for the three samples

10. Pages in your lab notebook containing the pertinent data

Questions on Cu in Brass Analysis

- 1. Why is it necessary to boil the water used to prepare the thiosulfate solution?
- 2. Why is Na₂CO₃ added to the thiosulfate solution?
- 3. Why is the thiosulfate solution stored in the dark?
- 4. Why is HCl added to the IO₃⁻ mixture and why must the solution be titrated immediately?
- 5. Why is the solution containing the dissolved brass sample heated to expel SO₃ fumes?
- 6. Why is H₃PO₄ added to the brass sample?
- 7. What is the purpose of the KSCN that is added just before the endpoint in the titration?

8. Why is the solution containing the dissolved brass made basic with concentrated NH_3 and then again acidified with H_2SO_4 ?

9. What is the formula of the tetrammine copper(II) complex?

10. Why do Zn^{2+} and Pb^{2+} not interfere in this procedure?

11. What sort of complications would arise if the iodinethiosulfate titration were carried out in a highly acidic solution?

12. If the solution were highly basic, how would the iodinethiosulfate reaction be influenced?

13. Why is the starch indicator not added at the beginning of the tritration?

Experiment 3: Estimation of Cu(II)

Aim:

Oxidation-Reduction (Redox) Reactions:

✤ Reaction involving change of oxidation number or transfer of electrons among the reacting substances.

The standard solutions are either oxidising or reducing agents.

The principal oxidising agents are KMnO4, K₂Cr₂O₇, Ce(SO₄)₂, I2, KIO₃, and KBrO₃.

♦ Frequently used reducing agents are Fe(II) and Sn(II) compounds, Na₂S₂O₃, As₂O₃, Hg₂(NO₃)₂, VCl₂ or VSO₄, CrCl₂ or CrSO₄ & TiCl₃ or Ti₂(SO₄)₃.

$2\mathbf{MnO}_{4}^{-} + 5\mathbf{H}_{5}\mathbf{C}_{5}\mathbf{O}_{4} + 6\mathbf{H}^{+}$	$\rightarrow 2Mn^{2+} + 10CO_2 + 8H_2$
(Pink)	(Colourless)
$2MnO_{4}^{-} + 10Fe^{2+} + 16H^{+}$	$\rightarrow 2Mn^{2+} + 10Fe^{3+} + 8H$
(Pink)	(Colourless)
$Cr_{2}O_{7}O_{7}O_{7}O_{7}O_{7}O_{7}O_{7}O_{7$	$\longrightarrow 2Cr^{3+} + 6Fe^{3+} + 7H_{,0}$
(Orange)	(Green)
CuSO ₄ + 4KI ——	\longrightarrow Cu,I, + K,SO ₄ +
(Blue)	(White ppt)
$2Na_2S_2O_3 + I_2$	$\longrightarrow Na_2S_4O_6 + 2NaI$
(Sod. Thiosulphate)	(Sod. Tetrathionate)

✤ The redox titration involving iodine directly or indirectly as an oxidizing agent are called iodine titrations.

Iodimetric Titrations: Titrations with a standard solution of iodine.

Applications of Iodimetric Titrations:

- ★ Determination of S2O3²⁻
- ★ Determination of SO3²⁻
- ★ Determination of AsO3³⁻
- ★ Determination of SnCl₂

 $2Na_2S_2O_3 + I_2 \longrightarrow Na_2S_4O_6 + NaI$ $Na_2SO_3 + I_2 + H_2O \longrightarrow Na_2SO_4 + 2HI$ $Na_3AsO_3 + I_2 + H_2O \longrightarrow Na_3AsO_4 + 2HI$

Iodometric Titrations: Iodine titrations in which some oxidizing agent liberate I_2 from an I⁻ solution and then liberated I_2 is titrated with a standard solution of reducing agent.

Estimation of CuSO_4 , $\operatorname{K}_2\operatorname{Cr}_2\operatorname{O}_7$, KMnO_4 , Fe^{3+} , $\operatorname{H}_2\operatorname{O}_2$, Br_2 , Cl_2 etc. $\begin{array}{c} \operatorname{CuSO}_4 + 4\operatorname{KI} & \longrightarrow & \operatorname{Cu}_2\operatorname{I}_2 + \operatorname{K}_2\operatorname{SO}_4 + \operatorname{I}_2 \\ (\operatorname{Blue}) & (\operatorname{White ppt}) \\ 2\operatorname{Na}_2\operatorname{S}_2\operatorname{O}_3 + \operatorname{I}_2 & \longrightarrow & \operatorname{Na}_2\operatorname{S}_4\operatorname{O}_6 + 2\operatorname{NaI} \\ (\operatorname{Sod. Thiosulphate}) & (\operatorname{Sod. Tetrathionate} \\ \operatorname{K}_2\operatorname{Cr}_2\operatorname{O}_7 + 4\operatorname{H}_2\operatorname{SO}_4 & \longrightarrow & \operatorname{Cr}_2(\operatorname{SO}_4)_3 + \operatorname{K}_2\operatorname{SO}_4 + 4\operatorname{H}_2\operatorname{O} + 3\operatorname{O} \\ \operatorname{Orange} & & \operatorname{Green} \\ 6\operatorname{KI} + 3\operatorname{H}_2\operatorname{SO}_4 + 3\operatorname{O} & \longrightarrow & 3\operatorname{K}_2\operatorname{SO}_4 + 3\operatorname{H}_2\operatorname{O} + 3\operatorname{I}_2 \\ 2\operatorname{Na}_2\operatorname{S}_2\operatorname{O}_3 + \operatorname{I}_2 & \longrightarrow & \operatorname{Na}_2\operatorname{S}_4\operatorname{O}_6 + 2\operatorname{NaI} \end{array}$

DETECTION OF THE END POINT:

- In this titration a solution of starch is used as indicator.
- Starch reacts with iodine in the presence of iodide to form an intensely blue-colored complex, which is visible at very low concentrations of iodine.
- The color sensitivity decreases with increasing temperature of the solution.
- It cannot be used in a strongly acid medium because hydrolysis of the starch occurs.
- The great merit of starch is that it is inexpensive.
- It is insoluble in cold water.
- Instability of suspensions in water
- It gives a water-insoluble complex with iodine, the formation of which precludes the addition of the indicator early in the titration (for this reason, in titrations of iodine, the starch solution should not be added until just prior to the end point when the colour begins to fade).
- There is sometimes a 'drift' end point, which is marked when the solutions are dilute.

Preparation and Use of Starch Solution:

- Make a paste of 0.1 g of soluble starch with a little water, and pour the paste, with constant stirring, into 100 mL of boiling water, and boil for 1 minute. Allow the solution to cool and add 2-3 g of potassium iodide. Keep the solution in a stoppered bottle.
- Only freshly prepared starch solution should be used.
- The same volume of starch solution should always be added in a titration.
- In the titration of iodine, starch must not be added until just before the end point is reached.
- Apart from the fact that the fading of the iodine colour is a good indication of the approach at the end point, if the starch solution is added when the iodine concentration is high, some iodine may remain adsorbed even at the end point.

The Starch-Iodine Complex

Numerous analytical procedures are based on redox titrations involving iodine.

Starch is the indicator of choice for these procedures because it forms an intense blue complex with iodine.

 \clubsuit Starch is not a redox indicator; it responds specifically to the presence of I₂, not to a change in redox potential.

♦ The active fraction of starch is amylose, a polymer of the sugar -D-glucose. Small molecules can fit into thecenter of the coiled, helical polymer.

 \clubsuit In the presence of starch, iodine forms I6 chains inside the amylose helix and the color turns dark blue.

I..I..I..I..I

Starch is readily biodegraded, so it should be freshly prepared.

 \diamond A hydrolysis product of starch is glucose, which is a reducing agent. Therefore, partially hydrolyzed starchsolution can be a source of error in a redox titration.

Preparation of Sodium thiosulphate Solution:

Sodium thiosulphate (Na2S2O3.5H2O) is readily obtainable in a state of high purity, but there is always some uncertainty as to the exact water content because of the efflorescent nature of the salt and for other reasons.

The substance is therefore unsuitable as a primary standard. It is a reducing agent by virtue of the half-cellreaction:

$$2S_2O_3^2 \longrightarrow S_4O_6^2$$

Recommendations for the Preparation of Sod. Thiosulphate:

◆ Prepare the solution with recently boiled distilled water.

Add 3 drops of chloroform, this improves the keeping qualities of the solution.

Avoid exposure to light, as this tends to hasten the decomposition.

Standardisation of Sod. Thiosulphate Solution:

 \diamond The standardisation of thiosulphate solutions may be done with potassium iodate, potassium dichromate, copper and iodine as primary standards, or with potassium permanganate or cerium(IV) sulphate as secondary standards.

Estimation of Cu(II) from Copper sulphate solution:

Pipette 10.0 mL of CuSO4 solution into a 250 mL conical flask then add a few drops of dilute sodium carbonate solution until a faint permanent precipitate remains, and this is removed by means of a drop or two of acetic acid. Then 0.25 g KI (or 2 mL of a 10 per cent solution) is added and titrate the liberated iodine with standard solution of sod. thiosulphate until the brown color of iodine fades, then add 0.4 mL of starch solution, and continue the addition of the thiosulphate solution until the blue colour commences to fade. Then add about 0.25 g of potassium thiocyanate or ammonium thiocyanate, preferably as a 10 per cent aqueous solution: the blue colour will instantly become more intense. Complete the titration as quickly as possible. The precipitate possesses a pale pink colour, and a distinct permanent end point is readily obtained.

Experiment

Aim: To estimate the strength in g/L of a given copper sulphate solution being provided with an approx. N/30 sod. thiosulphate (hypo) solution.

<u>Theory</u> : Step 1:	2CuSO	$_{4}^{+}$ 4KI \longrightarrow 2Cu ²⁺ + 2I ⁻	$Cu_{2}I_{2} + 2K_{2}SO_{4} + I_{2}$ $\longrightarrow 2Cu^{+} + I_{2}$
	01 .	2I ⁻ (aq.) 2Cu ²⁺ + 2e ⁻	$\xrightarrow{I_2 + 2e^-} 2Cu^+$
Step 2:	C	$I_2 + I^ $ Water Insoluble)	(Water soluble)
Step 3:		Starch + I ₃ (Colourless)	→ Starch-I ₃ (Starch-triiodide complex) (Dark Blue)
Step 4: Note:		$2S_2O_3^{2-} + I_2^{$	$\longrightarrow S_4O_6^{2-} + 2I^{-}$

Preparation of Standard solution:

 $W_1 = Weight of bottle + substance = _____ gms$

 $W_2 = Weight of bottle = _____ gms$

Weight of substance = $(W_1-W_2) = _____gms.$

Normality of the solution = $(W_1-W_2) \ge 10$ /Equivalent Weight =

 $= (W_1 - W_2) \times 10/49 =$

S.NO	Volume of the standard K ₂ Cr ₂ O ₇ solution (V1)	Burette Reading		
		final	Intial	Volume consumed
1	20 ml			
2	20 ml			
3	20 ml		2	

N1 = Normality of K2Cr2O7 solution =

V1 = Volume of K2Cr2O7 solution = 20ml

 $N_2 = Normality of Hypo =?$

 $V_2 = Volume of Hypo =$

 $N_1 V_1 = N_2 V_2$

 $N_2 = N_1 V_1 / V_2$ $N_2 = Normality of Hypo =$

S.No	Volume of the Copper solution (V3)	Burette Reading		Volume of consume (V4)
		Initial	final	
1	20ml			
2	20ml			
3	20ml			
				V4=

N3 = Normality of the Copper solution =? V3 = Volume of the Copper solution = 20ml N4 = Normality of Hypo =? V4 = Volume of Hypo = N3 V3= N4 V4 Normality of the Copper solution = N3 = N4 V4/V3

Amount of Copper present in the whole of the given solution $(100 \text{ ml}) = N_3 X 63.54/10$ **Result:** Amount of Copper present in the whole of the given solution $(100 \text{ ml}) = ______ \text{gm}$.

Why starch is added just before the end point?

The starch is not added at the beginning of the titration when the iodine concentration is high. Instead, it isadded just before the end point when the dilute iodine color becomes pale yellow. There are two reasons for such timing:First reason is that the iodine–starch complex is only slowly dissociated, and a diffuse endpoint would resultif a large amount of the iodine were adsorbed on the starch. The second reason is that most iodometric titrations are performed in strongly acid medium and the starch has a tendency to hydrolyze in acid solution.Starch-iodine complexation is temperature dependent. At 50 ^oC, the color is only one-tenth as intense as at 25 ^oC. If maximum sensitivity is required, cooling in ice water is recommended. Organic solvents decrease the affinity of iodine for starch and markedly reduce the utility of the indicator.

Why KSCN or NH4SCN is added?

When copper(II) is titrated iodometrically, the end point is diffuse unless thiocyanate ion is added. This isdue adsorption of I_2 on the surface of the cuprous iodide precipitate and only slowly reacts with thethiosulfate titrant. The thiocyanate coats the precipitate with CuSCN and displaces the iodine from the surface. The potassium thiocyanate should be added near the end point since it is slowly oxidized by iodineto sulfate. The pH must be buffered to around 3. If it is too high, copper(II) hydrolyzes and cupric hydroxide will precipitate. If it is too low, air oxidation of iodide becomes appreciable because it is catalyzed in the presence of copper.

Experiment 4: Estimation of Vitamin C

Vitamin C or ascorbic acid is essential for human life and is required for a range of physiological functions in human body. It can be found either in fresh fruits and vegetables naturally or in medical forms such as normal tablets, effervescent tablets and liquid vials. It is the most widely taken supplement. Though daily requirements of vitamin C are changeable according to the age, sex and conditions, it is around 75 to 90 mg per day for healthy adults and no more than 2000mg per day is recommended. It is one of the most ubiquitous vitamins ever discovered. Besides plays a paramount role as an antioxidant and free radical scavenger, it has been suggested to be an effective antiviral agent. In addition, ascorbic acid has been widely used in the pharmaceutical, chemical, cosmetic and food industry as antioxidant. Therefore, there is a need to find an accurate, reliable, rapid, and easy-to implement method for measuring the amount of ascorbic acid in a sample. However, there have been difficulties in quantifying ascorbic acid due to its instability in aqueous solution. The instability of ascorbic acid is due to its oxidation to dehydroascorbic acid, which is a reversible reaction, and subsequently to 2,3-diketo-L-gulonic acid. The later reaction is irreversible. Ascorbic acid is a water soluble vitamin with molecular weight of 176.12 g/mol and melting point of 193°C. World-wide accepted daily requirement of ascorbic acid is about 60-95 mg4. Ascorbic acid is a reducing agent which reverses the oxidation in aqueous solution. Increased amounts of free radicals trigger the condition called oxidative stress which is kept under control by antioxidants. If there are not enough antioxidants some stress related diseases including hypertension, atherosclerosis, chronic inflammatory diseases and diabetes might occur. The following Iodometric titration is performed and the amount of vitamin C was evaluated in given tablet using iodometric titration method. In this method the

reaction between iodine and starch suspension, will indicate the endpoint by producing the blueblack product. The tri-iodide ions are quickly converted into iodide ions when ascorbic acid is present. However, when all of the ascorbic acid is oxidized, the excess iodide will react with starch and will result in blue-black color.

Equipment and apparatus:

- > Volumetric flask 250.00 ml (\pm 0.01 ml)
- > Beaker 500.00 ml (\pm 0.01 ml), 500.00 ml (\pm 0.01 ml), 100 ml (\pm 0.01 ml)
- ➢ Burette 50.00 ml (± 0.01 ml)
- Staduated cylinder 50.00 ml (± 0.01 ml)
- Erlenmeyer flask 250.00 ml (± 0.01 ml), Erlenmeyer flask 250.00 ml
- > Potassium iodide (KI) $15.00g (\pm 0.001 g)$
- Iodine powder 5g
- ➢ Balance
- Starch powder 0.25g
- Distilled water
- ➢ Heater
- ➢ Glass rod
- ➢ Vitamin C

EXPERIMENTAL PROCEDURES:

Ascorbic acid is determined by using an oxidation-reduction reaction. The solubility of iodine is increased with iodide and tri-iodide is occurred:

$$I_2(aq) + I^- \leftrightarrow I_3^-$$

I₃ - then oxidizes vitamin C to dehydroascorbic acid:

 $C_6H_8O_6 + I_3 + H_2O \longrightarrow C_6H_6O_6 + 3I^- + 2H^+$ Vitamin C dehydroascorbic acid

The endpoint is production of a blue-black color which occurs as a result of the reaction of iodine with starch suspension. When ascorbic acid is present, I_3^- is converted to iodide and no color change is observed. However, when all ascorbic acid was utilized, expected blue-black color occurs due to the reaction between starch and excess tri-iodide.

This titration procedure is widely accepted and is appropriate for testing the amount of vitamin C in the tablets, liquids and fruits and vegetables.

Preparation of iodine solution:

For preparation of 0.1 M iodine solution, 10g of KI was taken in a 250 ml Volumetric flask and 35 ml of distilled water was added followed by heating the solution; the mixture was cooled to room temperature and 3.15 g of solid Iodine powder was dissolved.

Similarly, to prepare 0.005M of iodine solution 2g of KI was taken in a 500 ml beaker and dissolving in 100 ml of distilled water and 1.3 g of iodine powder was stirred with small quantity of water and qs (quantum satis) to 1 liter.

Preparation of starch solution:

Addition of 0.25 g of starch powder in 50ml warm distilled water, As the starch is insoluble in cold water and needs to be boiled to stay in solution.

Preparation of vitamin C standard solution:

25mg Ascorbic acid was taken in a 100.00 ml beaker and dissolved in 100 ml distilled water.

Preparation of Vitamin C sample solution:

From the strip of Vitamin C random two tablets were weighed and smashed to form powder and average value was calculated.

Weight of two tablets = 1.67g

Average = 1.67/2 = 835 mg

250 mg Ascorbic acid tablet has equivalent weight of = 835 mg

100mg Ascorbic acid tablet has equivalent weight of = 167*2 mg

= 334 mg

334 mg of the powder tablet was taken in a 100.00 ml beaker and dissolved in 100 ml distilled water.

Standardization of the iodine solution with the vitamin C standard solution and sample solution:

The measured volume of 20ml of both standard and sample was taken from each solution and equilibrated with 150 ml distilled water separately into distinct two Erlenmeyer flask 250.00 ml and titrant containing iodine solution was run against analyte containing either sample or standard; 5-6 drops of prepared starch solution were added to the analyte and titration was started.

The burette level for each analyte for distinctive experiment was noted as mentioned below: For standard solution the volume of iodine solution required for complete reaction = 45 ml Equally, for Sample solution the volume of iodine solution required = 49 ml The endpoint was noted when analyte appears blue in color.

CALCULATION

For sample solution

In the beginning of the experiment 20 ml of sample was taken from 100 ml of prepared solution containing 100 mg of Ascorbic acid. As 49 ml of iodine is required for the color change containing 20 ml ascorbic acid solution, the dilution was done 5 times to that of the solution. Hence, the final volume of the iodine solution = 49×5

 $Mole_{iodine} = Mass_{Ascorbic acid} \times \frac{1 \text{ mole}}{176.12g \text{ of ascorbic acid}} \times \frac{1000 \text{ml}}{\text{volume of iodine}}$

 $M_{iodine} = 0.1g \times 1/176.12 \times \frac{1000}{245 \, ml}$

For standard solution

Mass Ascorbic acid = Mole iodine \times Volume of iodine \times 176.12

 $= 0.00231 \times 45 \times 176.12$

= 91.54 mg

Initially, the amount of Ascorbic acid was taken for 100mg and therefore for total amount of ascorbic acid i.e. 250 mg the ratio stands out to be 2.5 (250/100)

RESULT

Therefore, a 250 mg tablet of ascorbic acid from the ACI LIMITED (Nutrivit® C) contain = $2.5 \times 91.54 = 228.85$ mg

Amount of Ascorbic acid in ACI LIMITED (Nutrivit® C) is 8.46% less than the claimed value. A relationship among the pharmaceutical companies listed can be illustrated with a bar chart. The default amount of Ascorbic acid for all the companies in a tablet dosage form is 250mg

Pharmaceutical Company	Brand Name	Amount of Ascorbic	Reduced amount
		acid (mg)	(%)
ACI LIMITED	Nutrivit® C	228.85	8.46
Opsonin	Vasco	248	0.8
ACME	Cecon	232.9	6.84
BEXIMCO	Ascobex	220	12%

Experiment 5: Estimation of Cr and Mn in Steel.

Objectives

After studying and performing this experiment you should be able to:

• Explain the principle underlying simultaneous spectrophotometric determination of species with overlapping spectrum,

• Carry out absorbance measurements for solutions on a spectrophotometer and draw its UV-VIS spectrum,

• Determine the wavelength of maximum absorption in the spectrum and compute the corresponding molar absorption coefficient,

• Calculate the concentrations of chromium and manganese in a mixture from the absorption measurements of the mixture at two different wavelengths,

• Apply and adapt the method of simultaneous determination for other similar combination of ions. **Principle:**

On dissolution of alloy steel containing chromium and manganese, we get these elements as their ions; Cr3+ and Mn2+ ions respectively. In the determination, these are oxidised respectively to dichromate and permanganate ions using potassium persulphate and potassium periodate, respectively. The orange red coloured dichromate shows maximum absorption (λ_{max}) at 440 nm while for the pink coloured permanganate the λ_{max} is at 545 nm. However, permanganate also absorbs at 440 nm to a smaller extent. Similarly, dichromate ions also have small absorption at 545 nm



Schematic diagrams of the visible spectra of potassium dichromate and potassium permanganate in 1M sulphuric acid, indicating the overlapping of the spectra and the relative values of molar absorption coefficients

In simple words, both the species absorb at the wavelengths of maximum absorptions mentioned above. Thus, the absorption of the solution containing a mixture of dichromate and permanganate ions at these wavelengths would be the sum of the absorbances of the two species. It can be shown that the absorbances of these two ions obey Beer-Lambert's law when measured individually and also in a mixture in presence of 0.5 M sulphuric acid. Under these conditions we can write the following expression assuming a unit path length.

$$A_{440} = \varepsilon_{Cr, 440} [Cr_2 O_7^{2-}] + \varepsilon_{Mn, 440} [Mn O_4^{-}] \qquad \dots (.1)$$

$$A_{545} = \varepsilon_{Cr, 545} [Cr_2 O_7^{2-}] + \varepsilon_{Mn, 545} [Mn O_4^{-}] \qquad \dots (.2)$$

Here, the meaning of the terms used is self evident. For example, ε Cr, 545 refers to the molar absorption coefficient at 545 nm for chromium as dichromate ion; A_{440} refers to the absorbance of the mixture at 440 nm and [Cr₂O₇²⁻] represents the concentration of dichromate ions and so on. We need to solve these simultaneous equations to determine both the species without their separation.

Solving Eq. (1) and (2) we get the following expressions for the concentrations of dichromate and permanganate ions. Thus, the concentrations of the two ions can be obtained from the absorbance measurements at 440 and 545 nm, if we know the molar absorption coefficients for both the ions at these wavelengths.

$$\begin{bmatrix} \operatorname{Cr}_{2}\operatorname{O}_{7}^{2-} \end{bmatrix} = \frac{A_{440}(\varepsilon_{\operatorname{MnO}_{4}})_{545} - A_{545}(\varepsilon_{\operatorname{MnO}_{4}})_{440}}{(\varepsilon_{\operatorname{Cr}_{2}\operatorname{O}_{7}^{2-}})_{440}(\varepsilon_{\operatorname{MnO}_{4}})_{545} - (\varepsilon_{\operatorname{MnO}_{4}})_{440}(\varepsilon_{\operatorname{Cr}_{2}\operatorname{O}_{7}^{2-}})_{545}} \qquad \dots (3)$$
$$\begin{bmatrix} \operatorname{MnO}_{4}^{-} \end{bmatrix} = \frac{A_{545} - (\varepsilon_{\operatorname{Cr}_{2}\operatorname{O}_{7}^{2-}})_{545}}{(\varepsilon_{\operatorname{MnO}_{7}})_{545}} \qquad \dots (4)$$

The values of the molar absorption coefficients are obtained from absorbance measurements of pure solutions of the two substances at the designated wavelengths. Fe^{3+} , Ni^{2+} , Co^{2+} and V^{2+} also absorb in these regions, therefore these should be absent from the analyte solution.

Requirements:

Apparatus	
Spectrophotometer/ Filter photometer	1
Matched cuvette	2
Volumetric flasks (1 litre)	1
Volumetric flasks (250 cm3)	2
Graduated pipette (10 cm3)	1
Beakers (100 cm3)	10
Weighing bottle	1
Burettes	2

Chemicals
Potassium dichromate
Potassium permanganate
Sulphuric acid

Solutions Provided:

i) **1 M sulphuric acid** (**Reference blank**): It is prepared by adding (with constant swirling) 53.3 cm3 of analytical grade concentrated sulphuric acid into a one litre flask containing about 500 cm3 of distilled water, allowing the flask to cool and making up the volume using distilled water.

ii) **Potassium dichromate (0.002 M) in 1M sulphuric acid**: It is prepared by accurately weighing 0.1471 g of analytical grade potassium dichromate and quantitatively transferring to a 250 cm3 volumetric flask and adding sufficient amount of 1 M sulphuric acid to dissolve it. Once the compound dissolves the volume is made upto the mark with 1M sulphuric acid.

iii) **Oxalic acid** (0.005 M): It is prepared by accurately weighing 0.063 g of oxalic acid and quantitatively transferring to a 100 cm3 volumetric flask and adding sufficient amount of water to dissolve it. Once dissolved the volume is made upto the mark with distilled water.

iv) Potassium permanganate (0.002 M) in 1 M sulphuric acid: It is prepared as follows.

- Weigh about 1.8 g A.R. potassium permanganate and dissolve in 50 cm3 of distilled water taken in a 100 cm3 beaker.
- Heat the solution to boiling and gently continue for about fifteen minutes.
- Cool the solution and filter through a plug of glass wool in the funnel.
- Take 10 cm3 of the filtrate in a 100 cm3 beaker and dilute to 100 cm3 with 1M sulphuric acid.
- Standardise the solution by titrating with a standard solution of oxalic acid in acidic medium.
- Calculate the volume of the solution required to prepare 250 cm3 of 0.002 M solution.
- Take the requited volume of the standardised solution of potassium permanganate in a 250 cm3 volumetric flask. Make up the volume with 1M sulphuric acid.

v) **Sample solution**: prepared by suitably mixing the 0.002 M solutions of potassium dichromate and potassium permanganate to make 10 cm3 of the mixture and adding 10 cm3 of 1M sulphuric acid to it.

Procedure:

The procedure for the experiment can be divided into three components. These are given below.

a) Determination of molar absorption coefficients of dichromate and permanganate ions

b) Establishing the additivity of absorbance values of dichromate and permanganate ions

c) Determination of the concentrations of dichromate and permanganate ions in the mixture.

Follow the instructions given below in the sequential order to accomplish these tasks.

a) Determination of molar absorption coefficients for dichromate and permanganate ions

1. Pipette out 20 and 40 cm3 respectively of the 0.002 M stock solution of potassium dichromate into two 100 cm3 standard flasks and make up the volume with 1M sulphuric acid. This would give 0.0004 M and 0.0008 M solutions respectively of potassium dichromate. Mark the flasks accordingly.

2. Pipette out 20 and 40 cm3 respectively of the 0.002 M stock solution of potassium permanganate into two 100 cm3 standard flasks and make up the volume with 1M sulphuric acid. This would give 0.0004 M and 0.0008 M solutions respectively of potassium permanganate. Mark the flasks accordingly.

3. Measure the absorbance values at 440 nm and 545 nm for the stock solutions (0.002 M) of potassium dichromate and potassium permanganate and the dilutions prepared in step 1 and 2. Make measurements in cuvettes of 1 cm path length and use 1M sulphuric acid as the reference. Record the values in Observation Table 1.

4. Calculate the molar absorption coefficients for potassium dichromate and potassium permanganate at 440 nm and 545 nm as per the expressions given under column 4 and 6 and record appropriately in the Observation Table 1.

b) Establishing the additivity of absorbance values of dichromate and permanganate ions at 440 and 545 nm

1. Take two clean burettes and fill them with 0.0008 M potassium dichromate and 0.0004 M potassium permanganate respectively.

2. Take 7 beakers of 100 cm3 capacity labelled from 1 to 7 and transfer the solutions of 0.0008 M potassium dichromate, 0.0004 M potassium permanganate and 1M sulphuric acid as per the details given in column 2, 3 and 4 of the Observation Table 2.

3. Measure the absorbance values of the resulting solutions using 1M sulphuric acid as the reference and record your observations in the column 5 and 7 of Observation Table 2. You must use cuvettes having a path length of 1 cm.

4. Calculate the absorbance values for the mixtures of potassium dichromate and potassium permanganate at 440 nm and 545 nm by using Eq. (1) and Eq. (2) respectively and record them in the column 6 and 8 of Observation Table 2.

5. Compare the observed and calculated absorbance values for the validation of additivity. Record your observations in the space provided after Observation Table 2.

c) Determination of the concentrations of dichromate and permanganate ions in the mixture.

1. Measure the absorbance values of the given sample solution taken in a cuvette of path length equal to 1 cm. Use 1M sulphuric acid as the reference and record your observations in the column 1 and 2 of Observation Table 3.

2. Calculate the concentration of potassium dichromate and potassium permanganate in the mixture by using Eq.(3) and Eq. (4) respectively and record your results in the column 3 and 4 of Observation Table 3.

Observations and Calculations:

A. Determination of molar absorption coefficients for dichromate and permanganate ions Observation Table-1: Determination of molar absorption coefficients for K₂Cr₂O₇ and KMnO4 at 440 and 545 nm

A: Potas	ssium dichromate	1.4.1	V les		
Column		2.4)			~
1	2	3	4	5	6
S. No.	Concentration of K ₂ Cr ₂ O ₇ c _{K₂Cr₂O₇}	Absorbanc e at 440 nm $= A_{440nm}$		Absorbance at 545 nm $= A_{545 nm}$	$\varepsilon_{545 \text{ mm}} \#$ $= \frac{A_{545 \text{ mm}}}{C_{\text{K}_2\text{Cr}_2\text{O}_7}}$
1	0.002				
2	0.0004				_
3	0.0008				
			Average =		Average =
B: Potas	sium permangana	te	50 		~
S. No.	Concentration of KMNO ₄ C _{KMnO4}	Absorbance at 440 nm $= A_{440nm}$	$\begin{split} \varepsilon_{440\mathrm{mm}} & \# \\ &= \frac{A_{440\mathrm{mm}}}{c_{\mathrm{KMnO_4}}} \end{split}$	Absorbance at 545 nm = A _{545 nm}	$\varepsilon_{545 \text{ mm}} \#$ $= \frac{A_{545 \text{ mm}}}{c_{\text{KMnO}_4}}$
1	0.002				
2	0.0004				
3	0.0008		ήr.		
-			Average =		Average =

The expression is valid only if the absorbance measure is made in a cuvette of path length=1cm.

The molar absorption coefficients for the $[Cr_2O_7^{2-}]$ and $[MnO_4^{-}]$ ions are found to be:

 $\epsilon_{440 \text{ nm}} [\text{Cr}_2 \text{O}_7^{2-}] =$ $\epsilon_{545 \text{ nm}} [\text{Cr}_2 \text{O}_7^{2-}] =$ $\epsilon_{440 \text{ nm}} [\text{Mn O}_4^{-}] =$ $\epsilon_{545 \text{ nm}} [\text{Mn O}_4^{-}] =$

B. Establishing the additivity of absorbance values of dichromate and permanganate ions at 440 and 545 $\rm nm$

Observation Table 2: The observed and calculated absorbance values for the mixtures of potassium dichromate and potassium permanganate.

Column							
1	2	3	4	5	6	7	8
S. No.	Volume of 0.0008 M K ₂ Cr ₂ O ₇	Volume of 0.0004 M KMnO4	Volume of 1 M H ₂ SO ₄	Absor at 440 $= A_{440}$	bance) nm	Abso at 54 $= A_{5}$	orbance 5 nm 45nm
	(cm [°])	(cm°)	(cm [°])	Obsd.	Calc.*	Obse	l. Calc.**
1.	50	0	1		1		
2.	40	10	1				
3.	30	20	1	9			- F.
4.	25	25	1				
5.	20	30	1				
6.	10	40	1				
7.	0	0	1				

* The absorbance value is calculated as per Eq. (3.1.)

** The absorbance value is calculated as per Eq. (3.2.)

Write your conclusions about the additivity of the absorbance values here.

.....

.....

C. Determination of the concentrations of dichromate and permanganate ions in the mixture. Observation Table 3: Determination of the concentrations of dichromate and permanganate ions in the mixture

Absorbance at 440 nm = A _{440nm}	Absorbance at 545 nm = A _{545nm}	$= \frac{C_{K_2Cr_2O_7}}{(\varepsilon_{Cr_2O_7^2})_{440} (\varepsilon_{MnO_4^2})_{545} - A_{545}(\varepsilon_{MnO_4^2})_{440}}$	$=\frac{C_{\text{KMnO}_4}}{(\varepsilon_{\text{MnO}_4})_{545}^2 - (\varepsilon_{\text{Cr}_2\text{O}_7^2})_{545}^2 \left[\text{Cr}_2\text{O}_7^{2-1}\right]}{(\varepsilon_{\text{MnO}_4})_{545}^2}$

 $C_{K_2Cr_2O_7} =$

 $C_{\text{KMnO}_4} =$

RESULT

i) The concentrations of potassium dichromate in the given mixture is found to be

=.....M

ii) The concentrations of potassium permanganate in the given mixture is found to be

=M

PRECAUTIONS

• All the reagents used must be of AR grade.

• All the glassware used must be quite clean.

 \cdot The measurements should preferably be made in cuvettes of 1 cm path length. (In the absence of these the expressions given for calculation of molar absorption coefficients and concentrations would need to be modified.)

Experiment 6: Estimation of (i) arsenite and (ii) antimony in tartar-emetic iodimetrically

Objectives:

1. The student will titrate using iodine as a titrant.

2. The student will use a buffer to control the reaction of the titration.

- 3. The student will use both acid and base to get the desired pH.
- 4. The standard and unknown are treated in the same manner.
- 5. Starch will be used as an indicator.

Preparation of Iodine Solution: Dissolve approximately 25 g of potassium iodide in about 50 mL of water. Add to this 6.5 g of iodine crystals. When the iodine has dissolved completely, dilute the solution to about500 mL. Do in a 600 mL beaker, not a volumetric flask!

Preparation of Starch Solution: Dissolve 1 gram of soluble starch in about 5 mL of cold d.i. Slowly add the starchsuspension to 95 mL of rapidly boiling water. Boil until the solution clears. Cool to roomtemperature. NOTE: Starch solutions often do not keep. Make it the day you will need it (ifyou also will need it the following period try to keep it--if it has a mold growth then it mustbe discarded).

Standardization of the Iodine Solution: Weigh out into 250 mL Erlenmeyer flasks FOUR samples of pure As2O3 of about 0.20 geach. Add 10 mL of 4% NaOH solution. Swirl until all As2O3 has dissolved. (You maywarm it to hasten solution, if necessary) Dilute to 50 mL. Cool to room temperature andadd 1:5 HCl until the solution turns pink with methyl red and cool again. Add SOLIDsodium bicarbonate in small portions until you are sure that no more will dissolve; thenadd about 3 g in EXCESS. If the bicarbonate is dissolved before a titration is completed,add more bicarbonate. Add 3 mL of starch solution and titrate with the iodine solutionimmediately. Continue the titration to the first permanent blue color. Calculate thenormality of the iodine solution.

Titration of Unknown Arsenic: Dry the unknown for 1 hour at 105 °C. Do not OVERHEAT. Cool it in a desiccator. Weighout into 250 mL Erlenmeyer flasks four samples of 0.4 to 0.5 g each. Add 10 mL of 4%NaOH solution. It may be necessary to warm it in order to hasten solution. Dilute to 50 mL.Cool to room temperature and add 1:5 HCl until solution turns pink with methyl red. Coolagain. Add bicarbonate in small portions until you are sure no more will dissolve; then addabout 3 g in excess. If the bicarbonate is dissolved before a titration is completed, add morebicarbonate. Add 3 mL of starch indicator and titrate IMMEDIATELY with the standardiodine solution. Calculate the percentage of As2O3 in the unknown. PRECAUTIONS:

1. Arsenic is a very poisonous element. Fatal dose is the amount that you can balance on thehead of a pin. Wash carefully before eating.

2. It is important that there be an excess of sodium bicarbonate throughout a titration.

3. Solutions of trivalent arsenic may be oxidized appreciable by air if allowed to stand. Thetitration should be completed as soon as possible after the sample is dissolved in alkali.

NOTES: Use 250 mL Erlenmeyer flasks for the titrations. Do not use aliquots.

The initial sample of arsenic oxide should be dried for one hour at 1000 C.

CALCULATION: Calculate the percent of arsenic oxide (As2O3) present in the sampleafter calculating the molarity of the titrant from the known arsenic titrations. Include aspread sheet and computer print-out for the calculation.

SAFETY AND DISPOSAL INFORMATION: Dispose of all solutions in the "As WASTE" container. Do not mix the wastes from the I2 and KMnO₄ Determinations. CAUTION: Remember to treat all chemicals with respect. Arsenic is BAD!

Experiment 7: Estimation of Fe in cement.

Aim: Estimation of Fe in cement colorimetrically. The method is based on the reaction of iron with ammonium thiocyanate forming a complex having a maximum absorption at 450 nm. Iron is estimated, the results are found to be reproducible and the method is economical and less time consuming. The method is successfully applied to cement sample and results is found to be accurate and precise as the more sophisticated colorimeter is commonly used for ferrous determination.

Composition of Cement

There are seven major ingredients of cement. The general percentage of these ingredients in cement is given below in Table 1:

Iron oxide: Chemical formula is Fe₂O₃.

 \Box Iron oxide imparts color to cement.

 \Box It acts as a flux.

 \Box At a very high temperature, it imparts into the chemical reaction with calcium and aluminum to form tricalcium alumino-ferrite.

□ Tricalcium alumino-ferrite imparts hardness and strength to cement

Colorometric measurements were made using a colorimeter and samples were analysed at 450 nm .The basis of spectrophotometric methods is the simple relationship between the absorption of radiation by the solution and the concentration of species in the solution. When monochromatic light passes through a transparent medium (coloured solution) the rate of decrease in intensity with the concentration and thickness of the medium is directly proportional to the intensity of the light .In order to determine a species or analyte in the solution spectrophotometrically, it is usually converted into a colored complex. Ammonium thiocyanate yields a blood red colour with ferrous iron and the colour produced is stable in nitric acid medium. The sample of cement was analysed with certain colouring reagent and optical density of the coloured compound is measured at 450 nm since the maximum absorption

was at 450nm, hence all measurements were made at 450 nm. Its optical density is measured in a colorimeter and the concentration of ferrous iron is found from the standard calibration curve. Methodology:

Dissolve the given ferrous ammonium sulphate (0.7022 gms) in 100 ml of water and add 5 ml of 1 : 5 H2SO4 followed by dil. KMnO4 solution through burette until light-pink colour appears. Dilute the solution to 1 L such that 1 ml of solution contains. 0.1 mg of Fe⁺². From the above

solution,	take	Ingredient	Percentage in cement	l
1,2,3,4,5	ml	Lime	60-65	into
100	ml	Silica	17-25	l
volumetri	с	Alumina	3-8	flasl
ml of	nitric	Magnesia	1-3	acid
of	40%	Iron oxide	0.5-6	
thiogyana	+070	Calcium Sulfate	0.1-0.5	ر برامی
the		Sulfur Trioxide	1-3	solu
uie a	anove	R.	•	sam

separately into separate standard flasks. Add 1 acid and 5 ml ammonium solution to all samples to get

blood red colour and make up the solutions to the mark by adding distilled water. Before recording the optical densities values for the above prepared solutions filter selection is done to know the maximum absorption in table 2.1. Now measure the optical densities of all the solutions using colorimeter table 2.2 and plot a graph by taking amount of ferrous iron on x-axis and optical density on the y-axis. The curve obtained is called standard calibration curve in figure 2.1

Dissolution of sample

Weight out 0.1 gm of the cement sample into a beaker. Add 5 ml of water and stir with a stirrer. Add few drops of concen HCl through the walls of the beaker and stir again. Heat the mixture until the moisture is evaporated. Then add 20 ml of distilled water and dissolve the content and make up the solution with distilled water to 100 ml. Shake well for uniform concentration.

Development of Colour

Pipette out 10 ml of above solution, into 100 ml std. flask add 1 ml of con HNO3 and then add 5 ml of 40% NH4SCN with a Burette. Make up the solution upto the mark with distilled water. Now take the sample solution into colorimeter tube and measure the O.D. (optical density) using the photo colorimeter. From the O.D. we can measure the concentration of Fe2+ from the calibration curve. The calibration curve is as follows in figure 1.

Filter Number	Absorbance
45	0.52
47	0.41
51	0.34
52	0.29
54	0.21
57	0.12

Table 1. Filter selection:

67	0.01

Table 2. Absorbance of solutions with va	rious concen of Ferrous(Fe ⁺²)
--	--

S.No	Concentration of Fe+2 in	Absorbance
	mg	
1	0.2	0.18
2	0.4	0.39
3	0.6	0.58
4	0.8	0.81
5	1	0.98



Figure 1 CALIBRATION CURVE

Calculation: From the calebration curve we determined the concentration of iron in cement.

C7P: Organic Chemistry Lab

Background

Organic chemists often must identify unknown compounds. In some cases, such as a reaction, you may have a good idea of what the compound in question is. However, in other cases, such as when you isolate a compound from a natural source, you may have no idea what the compound might be. In this experiment you will determine the identity of an unknown compound. First, you will need to purify your compound, then you will need to identify its functional group (it will contain only one), and finally you will need to make a derivative of the compound. You will confirm your results with boiling or melting point, IR, and NMR.

Impurities in your compound will make it extremely difficult to identify. Thus, before you do anything else, you will need to make sure your unknown compound is pure. Consider each of the following purification techniques you have learned over the course of the year.

- 1. *Recrystallization:* Works well for solid compounds. You will need to find an appropriate recrystallization solvent. Consider a variety of solvents and mixed solvent systems.
- 2. *Distillation:* Works well for liquids that have a boiling point of <250 °C. (Note: Fractional distillation may be required if you suspect impurities close to the boiling point of your unknown.)
- 3. *Column Chromatography:* Works well for UV active compounds. You will need to use TLC to identify a solvent system that will separate your unknown from any impurities.

After you have purified your unknown, verify that it is pure enough to proceed by measuring the boiling or melting point. Note that while you will not know what the melting point or boiling point of your unknown should be, the narrowness is an excellent indicator of whether or not your product is pure. Also pay attention to the appearance of your unknown and see if it has changed (hopefully for the better) during the course of the purification process.

Once your unknown is pure, you will need to identify its functional group. Your unknown will have one major functional group (alcohol, ketone, aldehyde, amide, amine, carboxylic acid, or ester). Additionally, your unknown compound may or may not contain an aromatic ring. To determine the functional group, it is recomanded that you start with solubility tests, and then conduct functional group classification tests. IR spectroscopy may also be useful at this point. Solubility can sometimes provide a surprisingly useful amount of information. First, you will test your unknown's solubility in water. Compounds with 4 carbons

or less will easily dissolve in water, whereas compounds with 8 carbons or more will be insoluble. Compounds containing 5----7 carbons may or may not dissolve (often they will display "partial" solubility). If your compound dissolves in water, you will also want to check the pH of the solution. Amines will typically be basic, and carboxylic acids will typically be acidic. Most other compounds will be neutral. Compounds that are insoluble in water should then be subjected to a solubility test in 5% HCl. Typically, only amines will be soluble in HCl because they form water-soluble hydrochloride salts when they react with HCl. Compounds that are not soluble in HCl, should be subjected to testing in basic solutions (5% NaOH and 5% NaHCO₃). Both strong and weak acids (Carboxylic acids and phenols) will be deprotonated by NaOH to form water-soluble alkoxides. Only strong acids like carboxylic acids will react with NaHCO₃. Compounds that are not soluble in base should then be reacted with a very strong acid, sulfuric acid (note that in the case of sulfuric acid, "solubility" is also indicated by any type of reaction such as heat, gas generation, or a color change). Compounds that cannot become protonated by sulfuric acid at all (i.e., alkanes, alkyl halides, and aromatic carbons) will still remain insoluble. These solubility tests are summarized in the flow charts below.



The results from the solubility tests can significantly help in determining which classification tests should then be performed, or at least narrow down the list. By no means do you need to conduct all classification tests. In fact, you should do your best to select only tests that will provide you with additional information about your unknown and/or confirm results. Also, make sure that your glassware is clean and dry so you do not get any false positive or false negative results. Keep in mind that a negative result for a classification test provides useful information, so be sure to keep track of negative results as well as positive results. Also, for each classification test that you perform, be sure to run a blank, and one or more controls. These will help you to determine if a reaction actually occurred. A blank includes everything but the unknown, and a control includes a compound for which the outcome is known in place of the unknown. Controls can be positive (a compound you know will react) or negative (a compound that you know will not react). The classification tests are summarized in the table below.

Functional group	Test	Test no	Notes
Elemental analysis	Lassaigne test	C-1	Test for nitrogen, sulphur,
			halogens
Amine	Basicity test	C-2	Test for aromatic amine
	Bleaching powder test	C-3	
	Dye test	C-4	
Nitro	Reduction test	C-5	Test for aromatic nitro
	Muliken and barker test	C-6	group
Amido	Nitrous acid test	C-7	Test for amido group
	Hydrolysis test	C-8	
Phenolic –OH	Ferric chloride test	C-9	Test for phenolic-OH
	Back dye test	C-10	
Carboxylic acid	Bicarbonate test	C-11	Test for carboxylic acid
	Esterification test	C-12	
Aldehyde	Benedict test	C-13	Test for aldehyde
	Tollens test	C-14	
Ketone	2,4-Dinitrophenoyl	C -15	Test for ketone
	Hydrazine test		

At this point, you should be able to use your boiling or melting point data combined with the results of your functional group data to develop a hypothesis as to what your unknown might be or at least narrow down the list to only a few candidates. Note tha due to the accuracy (or lack thereof) of our thermometers, your boiling or melting points may be up to 15 °C lower than the literature values.

Once your functional group has been determined, you will prepare a derivative of your unknown. To prepare a derivative, you will select a suitable reaction that converts your unknown into a different functional group for which the boiling or melting point is known. This is particularly useful because compounds that have similar boiling or melting points will oftenhave derivatives that differ significantly in terms of boiling or melting point. You should then be able to identify your unknown using this information.

Finally, you can confirm the identity of your product using IR and NMR. Note that these measurements can be taken at any time during the course of the lab after you purified your product. In fact, it is recommended that you conduct them sooner rather than later as they may provide valuable information as to the identity of your unknown (e.g., IR may reveal your functional group).

Lab Notebook Preparation A

Before coming to lab on the first day of this experiment, the following items must be in your lab notebook:

- 1. Title of experiment
- 2. Date the experiment is to be performed
- 3. Outline of your plan for determining the identity of your unknown
- 4. Hazards of and appropriate precautions for the safe handling of unknown organic compounds
- 5. References

Lab Notebook Preparation B

Before coming to lab on the day you plan to prepare a derivative, the following items must be in your lab notebook:

- 1. Title of experiment
- 2. Date the experiment is to be performed
- 3. List of possible unknowns
- 4. The chemical reaction(s) you are attempting (with skeletal structures...R groups are okay if you do not know the identity of your unknown yet)
- 5. For each reaction you are attempting, include a table with information about your starting materials. Include molecular weight, molar equivalents, and mmoles to be used. For solids include grams. For liquids, include grams, density, and volume. For solutions, include the concentration and volume. (Note: You will not be able to completely fill in the table if you do not know the identity of your unknown yet. If that is the case, list whatever data you can.)
- 6. Any relevant physical properties (i.e., melting points or boilingpoints of possible unknowns and their derivatives)
- 7. Hazards of and appropriate precautions for the specific reaction(s) you are conducting
- 8. References

Safety Notes

• Assume that all unknowns are flammable and harmful by inhalation, ingestion, and skin absorption. Do not inhale their vapors and avoid contact with eyes, skin and clothing.

Directions

- 1. Purify your unknown using distillation, recrystallization, or column chromatography. It is recommended that purify the entire unknown provided so that you have enough pure material for all of the tests.
- 2. Measure the boiling or melting point of your unknown to confirm its purity.
- 3. Confirm with your instructor that the boiling or melting point you obtained for your unknown is within 15 °C of the reported literature value before proceeding.
- 4. Test the solubility of your unknown in water. (If your unknown is a solid, crush it into a fine powder.)
 - a. Add approximately 30 mg of your unknown to a test tube or smallvial.
 - b. Add 1 mL of water and shake vigorously for approximately 30 seconds. If the unknown appears to be soluble, test the pH of the solution and then skip to step 9.
- 5. Test the solubility of your unknown in 5% HCl. (If your unknown is a solid, crush it into a fine powder.)
 - a. Add approximately 30 mg of your unknown to a test tube or small vial.
 - b. Add 1 mL of 5% HCl and shake vigorously for approximately 30 seconds. If the unknown appears to be soluble, skip to step 9.
- 6. Test the solubility of your unknown in 5% NaOH. (If your unknown is a solid, crush it into a fine powder.)
 - a. Add approximately 30 mg of your unknown to a test tube or small vial.
 - b. Add 1 mL of 5% NaOH and shake vigorously for approximately 30 seconds. If the unknown appears to be insoluble, skip to step8.
- 7. Test the solubility of your unknown in 5% NaHCO₃. (If your unknown is a solid, crush it into a fine powder.)
 - a. Add approximately 30 mg of your unknown to a test tube or small vial.
 - b. Add 1 mL of 5% NaHCO₃ and shake vigorously for approximately 30 seconds.
 - c. Note whether your unknown is soluble or insoluble and then skip to step 9.

- 8. Test the solubility of your unknown in concentrated H₂SO₄. (If your unknown is a solid, crush it into a fine powder.)
 - a. Add approximately 30 mg of your unknown to a test tube or small vial.
 - b. Add 1 mL of concentrated H₂SO₄ and shake vigorously for approximately 30 seconds.
 - c. Note whether your unknown is soluble or insoluble. (Any indication of a reaction such as heat, gas generation, or a color change also indicates solubility.)
- 9. Conduct classification tests as needed. See directions for specific tests below.
- 10. Confirm the identity of your functional group with your instructor before proceeding.
- 11. Prepare one or more derivatives of your unknown. See directions for specific derivatives below.
- 12. Measure the melting point of any derivatives.
- 13. Confirm with your instructor that the melting point you obtained for your derivative is within 15 °C of the reported literature value.

Classification Tests

C-1 Elemental Analysis

This reaction tests for the presence of nitrogens, sulphur and halogens.

Safety Notes: Sodium can cause serious burns and the sodium-lead alloy may react violently with some substances. Wear gloves, avoid contact, and keep the sodium-lead alloy away from other chemicals.

Recommended Controls: butylamine, acetamide, bromobenzene

- 1. In the fume hood or under a snorkel, place 0.25 g of 10% sodium-lead alloy in a small, dry test tube held vertically by a clamp.
- 2. Melt the alloy with a Bunsen burner flame and continue heating until the sodium vapor rises about 1 cm up the tube.
- 3. Using a Pasteur pipet, add 2 drops of the unknown (or 10 mg of a solid) directly onto the molten alloy so that it does not touch the sides of the tube.
- 4. Heat gently to start the reaction, remove the flame until the reaction subsides, then heat the tube strongly for a minute or two, keeping the bottom a dull red color.
- 5. Let the tube cool to room temperature.
- 6. Dropwise add 1.5 mL of water and heat gently for a minute or so until the excess sodium has decomposed and gas evolution ceases.

7. Filter the solution through a Pasteur pipet with a small plug of cotton, wash the cotton with 1 mL of water, and combine the wash water with the filtrate. (Use a rubber bulb to expel any liquid that adheres to the cotton.) The filtrate should be colorless or just slightly yellow. If it is darker, repeat the fusion with stronger heating or more of the alloy.

To test for nitrogen:

- 1. Put 5 drops of the sodium fusion solution into a small test tube.
- 2. While stirring, add enough solid sodium bicarbonate, to saturate it (a little excess solid should be present).
- 3. Add 1 drop of this solution to a test tube containing 10 drops of PNB reagent (*p*-nitrobenzaldehyde in dimethyl sulfoxide) and note any color change.

To test for halogens:

- 1. Acidify 10 drops of the sodium fusion solution with dilute nitric acid.
- 2. Boil it gently under the hood for a few minutes.
- 3. Add a drop or two of 0.3 M aqueous silver nitrate, and note the color and volume of any precipitate that forms. (If a voluminous precipitate forms, let the precipitate settle and then remove the solvent using a pipet.)
- 4. Add 2 mL of 3 M aqueous ammonia to the solid, shake vigorously, and note your observations.
- 5. To test further for bromine and iodine, acidify 1 mL of the original sodium fusion solution with 1 M sulfuric acid, boil for a few minutes, and add 0.5 mL of dichloromethane and a then a drop of freshly prepared chlorine water. Shake and look for a color in the dichloromethane layer. To test for sulphur:
 - 1. Put 1ml of sodium-extract into a test tube.
 - 2. Add 1 ml dil.NaOH solution followed by 2-3 drops of sodium nitroprusside.

Interpretation: In the PNB test, a purple color indicates the presence of nitrogen (green indicates sulfur). In the halogen tests, formation of a voluminous precipitate on addition of silver nitrate indicates that a halogen is present, and the color of the precipitate (a silver halide) may suggest which halogen: white for chlorine, pale yellow for bromine, and yellow for iodine. If only a faint turbidity is produced, it may be caused by traces of impurities or by incomplete sodium fusion. If the precipitate is silver chloride, it will dissolve in aqueous ammonia; silver bromide is only slightly soluble and silver iodide is insoluble. In the chlorine water test, a red-brown color is produced by elemental bromine and a violet color by elemental iodine. In the sulphur test, a violet or purple color indicates the presence of sulphur.

C-2 Basicity Test

This test is useful if you have already determined that you have an amine. It is used to distinguish alkyl amines from aromatic amines.

Recommended Controls: p-toluidine, dibutylamine

Procedure for water-soluble compounds:

- 1. Dissolve 4 drops of your unknown (0.10 g of a solid) in 3 mL of water.
- 2. Measure the pH of the solution using pH paper.

Procedure for water-insoluble compounds:

- 1. Dissolve 4 drops of your unknown (0.10 g of a solid) in 3 mL of a pH 5.5 acetate-acetic acid buffer.
- 2. Mix thoroughly.

Interpretation: Water-soluble alkyl amines give pH values above 11, whereas water-soluble aromatic amines have pH values below 10. Water-insoluble alkyl amines should dissolve in the buffer, but water-insoluble aromatic amines will not dissolve.

C-3 Bleaching powder test

Procedure:

- 1. Dissolve 0.05 g of your unknown in 5 mL of water.
- 2. Add 3-4 drops of bleaching powder solution.
- 3. Shake vigorously.

Interpretation: Transient purple color which soon turns brown or light purple color.

C-4 Dye test:



Procedure:

- 1. 0.1 g of organic sample is dissolved in 5ml of dil.HCl.
- 2. The mixture is cooled at 0°-5°C in an ice-bath.
- 3. Then add 1ml of ice cold solution of dil. NaNO₂.
- 4. The mixture is added to ice-cold alkaline solution of β -Napthol.

Interpretation: Red or orange red dye (brown or raddish purple or violet dye indicates the presence of two amino groups; soluble dye indicates the presence of SO3H or Ar-OH along with Ar-NH₂ group).

C-5 Reduction test



- 1. A mixture of 0.1 g of organic sample, few pieces of granulated tin or zinc and 3ml of Conc. HCl is warmed gently with occasional shaking till the reaction is complete.
- 2. The mixture is cooled.
- 3. Filtered, if required, diluted and diazo-coupling reaction is performed.

Interpretation: Brilliant red or scarlet dye obtained.

C-6 Muliken and Barker Test

Procedure:

- 1. 0.1 g of organic sample is dissolved in 5 ml of 50% alcohol.
- 2. A little solid NH₄Cl or 10% CaCl₂ solution and a pinch of Zn-dust is added to it.
- 3. The mixture is boiled for a few minutes.
- 4. Then the mixture is cooled and allowed to stand for 5 minutes and then filtered.
- 5. With the filtrate following three tests are performed:
- a) A portion of the solution is added to Tollen's reagent and then warmed in a water bath.
- b) Two drops of benzoyl chloride and 2 drops of conc. HCl are added to another portion of the filtrate followed by 12 drops of FeCl₃ solution.
- c) The last portion of the filtrate is warmed with a little Fehling's solution.

Interpretation: From the part (a), a silver mirror or black or grey precipitation is obtained. From part (b), a wine-red color of ferric hydroxamate is present, from last part (c), a red precipitation is obtained.

C-7 Nitrous Acid Test

Procedure:

1. A little of the aqueous solution of organic sample is treated with with a few drops of HNO₂ (NaNO2 and HCl).

Interpretation: Effervescence due to evolution of N₂ gas.

C-8 Hydrolysis Test

Procedure:

0.2 g of organic sample is heated with 2ml of 50% NaOH solution.

Interpretation: Characteristics smell of NH₃ which turns mercurous nitrate paper black or copper sulphate paper deep blue.

C-9 Ferric Chloride

This reaction tests for the presence of phenols.

Recommended Control: phenol

Procedure:

- 1. Dissolve 1 drop of the unknown (40 mg of a solid) in 1 mL of water. If(you know based on the results of your solubility tests that the unknown is insoluble in water, use 0.5 mL of water and 0.5 mL of methanol instead of 1 mL of water.)
- 2. Add two drops of 2.5% ferric chloride solution.

Interpretation: Formation of an intense red, green, blue, or purple color suggests a phenol or an easily enolizable compound (such as an aldehyde or ketone). Some phenols do not react under these conditions.

C-10 Back Dye Test

This reaction tests for the presence of phenols.

Procedure:

- 1. A few drops of aniline dissolved in dil. HCl.
- 2. Few drops of cold dil. NaNO₂ solution is added.
- 3. Then the clear solution is added to the cold solution of organic sample in NaOH.

Interpretation: A brilliant red dye is obtained. Phenolic OH group present and confirmed.

C-11 Bicarbonate Test

This reaction tests for the presence of carboxylic acid.

Procedure:

1. A small amount of organic sample is sprinkled over aqueous solution of sodium bicarbonate.

Interpretation: Effervescence due to the evolution of CO₂.

C-12 Esterification Test

This reaction tests for the presence of carboxylic acid.

- 1. 0.5 g of organic sample is taken in a dry test tube.
- 2. To this, add 1 ml of dehydrated ethanol.
- 3. Then 2-3 drops of conc. H_2SO_4 is added and heated for 5 minutes in a water bath.
- 4. The mixture is then poured into a beaker containing large volume of Na₂CO₃ solution.

Interpretation: Characteristics sweet fruity smell of ester.

C-13 Benedict's Test

This reaction tests for the presence of aldehydes. Note that most ketones and aromatic aldehydes will not react.



Recommended Controls: butanal

Procedure:

- 1. Add 2 drops of the unknown (80 mg if it is a solid) to 2 mL of water.
- 2. Add 2 mL of Benedict's reagent.
- 3. Heat the mixture to a boil.
- 4. Observe if a precipitate forms, and note its color.

Interpretation: Benedict's reagent contains copper(II) sulfate, sodium citrate, and sodium carbonate. Aldehydes will react with the Cu_{2+} from the copper(II) sulfate to form copper(I) oxide which appears as a yellow or orange precipitate (it may look a little green in the blue reaction solution). Note that most ketones and aromatic aldehydes will not react.

C-14 Tollen's Test

This reaction tests for the presence of aldehydes.

Recommended Controls: benzaldehyde

Procedure:

- 1. Measure 2 mL of 0.3 M aqueous silver nitrate into a test tube and add 1 drop of 3 M sodium hydroxide.
- 2. Add 2 M aqueous ammonia drop by drop, with shaking, until the precipitate of silver oxide just dissolves (avoid an excess of ammonia).
- 3. Add 1 drop of the unknown (40 mg of a solid) to this solution, shake the mixture, and let it stand for 10 minutes. (If a silver mirror is observed at this point, this is considered a positive result.)
- 4. Heat the mixture in a 35 °C water bath for 5 minutes.
- 5. Immediately after the test has been completed, dissolve any solid residue in 1M nitric acid and then dispose of the solution in the designated waste container.
- 6. *Interpretation:* Formation of a silver mirror on the inside of the test tube is a positive test for an aldehyde. (Note that if the tube is not sufficiently clean, a black precipitate or a suspension of metallic silver may form instead.)

C-15 2,4-Dinitrophenylhydrazine

This reaction tests for the presence of aldehydes and ketones.

Safety Notes: 2,4-Dinitrophenylhydrazine (DNPH) is harmful if absorbed through the skin. Wear gloves and avoid contact.

Recommended Controls: cyclohexanone, benzaldehyde

Procedure:

- 1. Dissolve 1 drop of the unknown (40 mg of a solid) in 1 mL of 95% ethanol (use more ethanol if necessary to completely dissolve the unknown).
- 2. Add this solution to 2 mL of the DNPH reagent.
- 3. Shake and let the mixture stand for 15 minutes or until a precipitate forms. (If a precipitate is observed at this point, this is considered a positive result.)
- 4. Scratch the inside of the test tube and observe if a precipitate forms, and note its color.

Interpretation: Formation of a crystalline yellow or orange-red precipitate indicates an aldehyde or ketone. The color of the precipitate may give a clue to the structure of the carbonyl compound (unconjugated aliphatic aldehydes and ketones usually yield a yellow precipitate, while aromatic and α , β -unsaturated aldehydes and ketones yield a orange-red precipitate).

Derivatives of Alcohols

D-1 p-Nitrobenzoates and 3,5-Dinitrobenzoates

For these derivatives, it is extremely important to ensure that your glassware and your alcohol are dry (i.e., free of water). Water can easily react with the acid chlorides to form carboxylic acids to form the respective carboxylic acids rather than the desired esters. (Note: The *p*----nitrobenzoic acid has a melting point of 237 °C. and the 3,5----dinitrobenzoic acid has a melting point of 205----207 °C.) If necessary, dry your glassware in an oven before proceeding.



- 1. Dry your unknown alcohol with magnesium sulfate or sodium sulfate.
- 2. Filter to remove the drying agent.
- 3. If you are making the *p*-nitrobenzoate derivative, add 0.20 g of *p*-nitrobenzoyl chloride to a small round bottom flask. If you are making the 3,5-dinitrobenzoate derivative, add 0.20 g of 3,5-dinitrobenzoyl chloride to a small round bottom flask.
- 4. In the fume hood or under a snorkel, dropwise add 0.10 g of your unknown alcohol to the acid chloride while stirring.
- 5. Heat the mixture in a 60-70 °C water bath. If your alcohol has boiling point of < 160 °C, heat the mixture for 5 minutes. If your alcohol has boiling point of > 160 °C, heat the mixture for 15 minutes.
- 6. Stir in 4 mL of 0.2 M sodium carbonate.
- 7. Heat the mixture to 50-60 °C for 30 seconds.
- 8. Cool to room temperature and then in an ice bath.
- 9. Collect the precipitate by small-scale vacuum filtration. Wash with cold water.

10. Recrystallize the precipitate from ethanol or an ethanol-water mixture.

D-2 Phenylurethanes and 1-Naphthylurethanes

For these derivatives, it is extremely important to ensure that your glassware and your alcohol are dry (i.e., free of water). Water can easily react with the isocyanates to form the respective ureas rather than the desired carbamates. (Note: The dipheylurea has a melting point of 241 °C, and the di-1-naphthylurea has a melting point of 297 °C.) If necessary, dry your glassware in an oven before proceeding.



Safety Notes: Isocyanates are irritants and lachrymators. Avoid contact with these reagents and use in a fume hood or under a snorkel.

- 1. Dry your unknown alcohol with magnesium sulfate or sodium sulfate.
- 2. Filter to remove the drying agent.
- 3. If you are making the phenylurethane derivative, add 5 drops of phenyl isocyanate to a small round bottom flask. If you are making the 1-naphthylurethane derivative, add 5 drops of 1-naphthyl isocyanate to small round bottom flask.
- 4. In the fume hood or under a snorkel, dropwise add 5 drops of the dry alcohol to the isocyanate. If no reaction is apparent, heat the mixture in a 60-70 °C water bath for 15 minutes.
- 5. Cool to room temperature and then in an ice bath.
- 6. Collect the precipitate by small-scale vacuum filtration.
- 7. Recrystallize the precipitate from petroleum ether or heptane. (If necessary, perform a hot gravity filtration.)

Derivatives of Aldehydes and Ketones

D-3 2,4-Dinitrophenylhydrazones



Safety Notes: 2,4-Dinitrophenylhydrazine is toxic and sulfuric acid is corrosive. Avoid contact with these reagents and use in a fume hood or under a snorkel.

Procedure:

- 1. In a small round bottom flask, dissolve 0.10 g of the unknown aldehyde or ketone in 1 mL of ethanol. (If your unknown is not completely dissolved, add ethanol drop bydrop until it goes into solution).
- 2. In the fume hood or under a snorkel, dropwise add 3 mL of the 2,4-dinitrophenylhydrazinesulfuric acid reagent.
- 3. Allow the solution to stand at room temperature until crystallization is complete. If no crystals form, heat the mixture in a 60-70 °C water bath for 15 minutes, and then it cool again. If there is still no precipitate, add cold water drop by drop to the solution until precipitate is observed.
- 4. Collect the precipitate by small-scale vacuum filtration. Wash once with 5 mL of cold 5% NaHCO₃ and once with cold water.
- 5. Recrystallize the precipitate from ethanol or an ethanol-water mixture. (Note: If more than 6 mL of ethanol is needed for recrystallization, add ethyl acetate drop by drop to the hot solution until everything is dissolved.)

D-4 Semicarbazones



Safety Notes: Semicarbazide hydrochloride is a suspected carcinogen. Avoid contact with the reagent.

- 1. Mix together 0.20 g of semicarbazide hydrochloride, 0.30 g of sodium acetate, and 2 mL of water in a small round bottom flask.
- 2. If your unknown aldehyde or ketone is water soluble, add 0.20 g of it directly to the flask and stir to dissolve. If your unknown aldehyde or ketone is not water soluble, add a minimum amount of ethanol to the mixture until your unknown goes into solution.

- 3. Stir the mixture for two minutes.
- 4. Cool the mixture in an ice bath. If no crystals form, heat the mixture in a 60-70 °C water bath for 5 minutes, and then it cool again.
- 5. Collect the precipitate by small-scale vacuum filtration. Wash with cold water.
- 6. Recrystallize the precipitate from ethanol or an ethanol-water mixture.

D-5 Oximes



Safety Notes: Hydroxylamine hydrochloride is toxic and mutagenic. Avoid contact with the reagent.

Procedure:

- 1. Mix together 0.125 g of hydroxylaminehydrochloride, 0.30 g of sodium acetate, and 2 mL of water to a small round bottom flask.
- 2. If your unknown aldehyde or ketone is water soluble, add 0.20 g of it directly to the flask and stir to dissolve. If your unknown aldehyde or ketone is not water soluble, add a minimum amount of ethanol to the mixture until your unknown goes into solution.
- 3. Stir the mixture for two minutes.
- 4. Cool the mixture in an ice bath. If no crystals form, heat the mixture in a 60-70 °C water bath for 15 minutes, and then it cool again. If there is still no precipitate, add cold water drop by drop to the solution until precipitate is observed.
- 5. Collect the precipitate by small-scale vacuum filtration. Wash with cold water.
- 6. Recrystallize the precipitate from ethanol or an ethanol-water mixture.

Reporting Your Results

Write your report according to the guidelines described in "Topic 4: Writing an Organic Chemistry Lab Report". Work by yourself on this report.

References & Additional Resources

1. Lehman, J. W. *Operational Organic Chemistry: A Problem----Solving Approach to the Laboratory Course*, 3rd ed.; Prentice Hall: Upper Saddle River, NJ, 1999; pp 529----572.

Alcohol	BP (°C)	3,5-Dinitro benzoate MP (°C)	4-Nitro benzoate MP (°C)	1-Naphthyl urethane MP (°C)	Phenyl urethane MP (°C)
methanol	65	108	96	124	47
ethanol	78	93	57	79	52
2-propanol	83	123	110	106	88
2-methyl-2-propanol	83	142	-	-	136
2-propen-1-ol	97	49	28	108	70
1-propanol	97	74	35	105	57
2-butanol	99	76	26	97	65
2-methyl-2-butanol	101	116	85	72	42
2-methyl-1-propanol	108	87	69	104	86
3-pentanol	116	101	17	95	48
1-butanol	118	64	36	71	61
2-pentanol	119	62	24	74	-
3-methyl-3-pentanol	123	94 (62)	69	104	43
3-methyl-1-butanol	132	61	21	68	57
4-methyl-2-pentanol	132	65	26	88	143
1-pentanol	137	46	11	68	46
cyclopentanol	141	115	62	118	132
2-ethyl-1-butanol	148	51	-	60	-
1-hexanol	157	58	5	59	42
cyclohexanol	161	113	50	129	82
furfuryl alcohol	172	80	76	130	45
1-heptanol	177	47	10	62	60
2-octanol	174	32	28	63	114
1-octanol	195	61	12	67	74
1-phenylethanol	202	95	43	106	92
benzyl alcohol	204	113	85	134	77
2-phenylethanol	219	108	62	119	78
1-decanol	231	57	30	73	60
3-phenyl-1-propanol	236	45	47	-	92
1-dodecanol	259	60	45 (42)	80	74

Table 1. Possible Alcohol Unknowns

*A dash indicates that no information is reported in the literature

**Melting points in parenthesis represent conflicting literature values

Alcohol	BP (°C)	МР (°С)	2,4-Dinitro phenylhydrazone MP (°C)	Semicarbazone MP (°C)	Oxime MP (°C)
ethanal	21		168 (157)	162	47
propanal	48		148 (155)	154	40
propenal	52		165	171	-
2-methylpropanal	64		187 (183)	125 (119)	-
butanal	75		123	106	-
3-methylbutanal	92		123	107	48
pentanal	103		106 (98)	-	52
2-butenal	104		190	199	119
2-ethylbutanal	117		95 (30)	99	-
hexanal	130		104	106	51
heptanal	153		106	109	57
2-furaldehyde	162		212 (230)	202	91
2-ethylhexanal	163		114 (120)	254d	-
octanal	171		106	101	60
benzaldehyde	178		239	222	35
phenylethanal	195	33	121 (110)	153 (156)	99
4-methylbenzaldehyde	204		234	234 (215)	80
3,7-dimethyl-6-octenal	207		77	84 (91)	-
2-chlorobenzaldehyde	209		213 (209)	229 (146)	76 (101)
4-methoxybenzaldehyde	248		253d	210	133
2-methoxybenzaldehyde		38	254	215	92
4-chlorobenzaldehyde		48	265	230	110 (146)
3-nitrobenzaldehyde		58	290	246	120
4-nitrobenzaldehyde		106	320	221 (211)	133 (182)

*A dash indicates that no information is reported in the literature **Melting points in parenthesis represent conflicting literature values

***If the substance changes color and smokes, this is considered decomposition (d = decomposes)

Table 3. Possible Ketone Unknowns

Ketone	BP (°C)	MP (°C)	2,4-Dinitro phenylhydrazone MP (°C)	Semicarbazone MP (°C)	Oxime MP (°C)
acetone	56		126	187	59
2-butanone	80		118	146	-
3-methyl-2-butanone	94		124	113	-
2-pentanone	101		143	112 (106)	58
3-pentanone	102		156	138	69
3,3-dimethyl-2-butanone	106		125	157	75 (79)
4-methyl-2-pentanone	117		95 (81)	132	58
2,4-dimethyl-3-pentanone	124		95 (88)	160	34
2-hexanone	128		110	125	49
4-methyl-3-penten-2-one	130		205	164 (133)	48
cyclopentanone	131		146	210 (203)	56
4-heptanone	144		75	132	-
2-heptanone	151		89	123	-
cyclohexanone	156		162	166	91
2,6-dimethyl-4-heptanone	168		92	122	210
2-octanone	173		58	124	-
cycloheptanone	181		148	163	23
acetophenone	202		238	198 (203)	60
2-methylacetophenone	214		159	205	61
propiophenone	218	21	191	182 (174)	54
3-methylacetophenone	220		207	203	57
4-methylacetophenone	226	28	258	205	88
2-undecanone	228		63	122	44
4-phenyl-2-butanone	235		127	142	87
3-methoxyacetophenone	240		-	196	-
2-methoxyacetophenone	245		-	183	83 (96)
4-methoxyacetophenone		38	228	198	87
4-phenyl-3-buten-2-one		42	227 (223)	187	117
benzophenone		48	238	167	144
2-acetonaphthone		54	262d	235	145
3-nitroacetophenone		80	228	257	132
9-fluorenone		83	283	234	195

*A dash indicates that no information is reported in the literature

**Melting points in parenthesis represent conflicting literature values

***If the substance changes color and smokes, this is considered decomposition (d = decomposes)

GE-3: Lab

Identification of a pure organic compound

Liquid Compounds

1. Methyl Alcohol CH₃-OH

Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat 7. Ignition test

- : Liquid
- : Colorless
- : Pungent but faint alcoholic
- : Miscible With water
- : Neutral
- : Volatilises
- : Blue non-sooty flame

Special Tests

Experiment	Observation
1. Oxidation test : A Cu-spiral is made	1.
repeatedly red hot and introduced into 2 ml	
of O.S. kept in a test tube dipped in a beaker	
containing cold Water. The solution is then	
divided into four pan ts.	
(a) Schiffis test : One part is added to	(a) Pink color which deepens slowly.
Schiff s.	
(b) Tollen's test: Another part is added to	(b) Bright mirror of silver on the side of the
Tollen's reagent and heated in a water-bath.	test tube.
(c) Resorcinol test: A mixture of 0.5 ml of	
oxidised liquid and a drop of 0.5% aq.	(c) A reddish -violet ring appears at the
Solution of resorcinol is added carefully	junction of two liquid layers. A white ppt.

down the side of the test tube held in an	changing to reddish-violet appear in the aq.
inclined position containing 2ml of conc.	Layer above the ring after a little while.
H_2SO_4 .	
2. Denige's reagent: To about 5ml of aq.	2. A violet color develops deepening on
Solution of o.s. in a porcelain basin placed	standing.
on ice-water, 2-3 ml of 2.5% soln. of	
KMnO ₄ soln. is added followed by few	
drops of conc. H ₂ SO ₄ when a brown color	
developes. The brown color is destroyed by	
addition of saturated soln. of oxalic acid. A	
freshly prepared soln.of schiffs reagent is	
added dropwise to it with stirring.	
	3. Characteristics smell of methyl
3. Oil of wintergreen test: 0.5 ml of o.s.	salicylate.
and 3 drops of conc. H2SO4 are added to	
0.5 g of salicylic acid and heated for one	
minute and finally the reaction mixture is	
poured into 50 ml of water taken in a	4. The test is negative i.e., no yellowish
beaker.	ppt. of iodoform
4. Iodoform test in NaOH is performed	(distinction from ethyl alcohol and
(vide test 3a for ethanol)	acetone)

2. Ethanol CH₃CH₂OH

Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Liquid
- : Colorless
- : Pungent but faint alcoholic
- : Miscible with water
- : Neutral
- : Volatilises
- : Blue non-sooty flame

Special Tests

Experiment	Observation

1. Oxidation test : A Cu-spiral is made	1.
repeatedly red hot and introduced into 2 ml	
of O.S. kept in a test tube dipped in a beaker	
containing cold Water. The solution is then	
divided into four parts.	
(a) Schiff's test : One part is added to	(a) Pink color which deepens slowly.
Schiff s.	
(b) Tollen's test: Another part is added to	(b) A silver mirror is disposited on the
Tollen's reagent and heated in a water-bath.	walls of the test tube.
(c) To another part equal volume of very	
dilute solution of sodium nitroprusside	(c) Wine-red color develops.
solution is added followed by a few drops	
of NaOH solution.	
(d) To another part a very dilute solution of	
sodium nitroprusside solution is added	(d) A deep red color develops.
followed by a few drops of piperidine.	
2. Ethyl acetate formation: About 0.2 ml	
of O.S. and an equal volume of glacial	
acetic acid and a little conc. H_2SO_4 (or 0.2	
g of fused sodium acetate) are taken in a	2. Pleasant fruity smell.
test tube and warmed in a water-bath for	
about 5 min. the mixture is then poured into	
water taken in a beaker to which a little	
Na ₂ CO ₃ has been added.	
3. Iodoform test: (a) To 2 ml of aq. soln.	
of o.s. an equal volume of a conc. solution	
of iodine in potassium iodide is added and	2 (a) Vallow ppt of indeform and its
then dil. NaOH solution is added dropwise	s. (a) Tenow ppt. of fodotorin and its
with stirring until violet color of iodine	(distinction from moth an of)
disappears. Then the mixture is warmed	(aisunction from methanol)
and cooled under tap with shaking.	
(b) The iodoform test is repeated by adding	
NH4OH instead of NaOH.	(b) No ppt. of iodoform.
	(distinction from acetone)

3. Chloroform CHCl₃

Physical characteristics and preliminary test :

- 1. State 2. Colour
- 3. Odour

- : Liquid
- : Colorless
- : Sweet

- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Immiscible with water and much heavier than water : Neutral
- : Volatilises
- : Yellow sooty flame

Special Tests

Experiment	Observation
1. Hydrolysis test: O.S. is treated with aq. KOH or NaOH solution and boiled. The mixture is acidified with conc. HNO ₃ and then AgNO ₃ is added.	 Curdy white ppt. soluble in NH₄OH but reappears on addition of conc. HNO_{3.} A vellowish red ppt.
 2. Fehlings test: O.S and Fehlings solution (I + II equal vol.) is warmed with constant shaking. 3. Resorcinol test: A little powdered resorcinol and few drops of O.S. is taken in a basin, about 1 ml of conc. solution of NaOH is added to it, then the mixture is warmed gently. 	 3. A brilliant reddish colouration is developed in aq. layer. 4. Intolerable obnoxious smell of carbylamine.
4. Carbylamine test: Few drops of O.S. aniline and alc. KOH are warmed in a diy test tube and the ensuing gas is smelt by placing the thumb at the mouth of the test tube and then holding the thumb near the nostrils.	(Not recommended to perform)





Physical characteristics and preliminary test :

- 1. State
- 2. Colour
- 3. Odour
- 4. Miscibility
- 5. Litmus

- : Liquid
- : Pale yellow
- : Characteristic smell of bitter almonds
- : Immiscible with water
- : Neutral
6. Action of heat

7. Ignition test

: Volatilises without leaving any

- residue
- : Yellow sooty flame

Special Tests

Experiment	Observation
1. Caustic soda test : About 0.2 ml of	1. The color darkens.
O.S. is heated with conc. solution of	
caustic soda.	
2. Reduction and diazocoupling test: 0.5	
ml is reduced by tin or zinc and dil. HCl	
for 5 minutes. The decanted solution is	2. Brilliant scarlet-red dye.
cooled and very dilute cold solution of	
NaNO ₂ is added to it. The diazotised	
soln. is added to cold alkaline soln. of. β -	
naphthol.	3 Grev or black ppt
3. Muliken-Barker test: An aq. ethanolic	5. Grey of black ppt
solution of few drops of given sample is	
boiled with a pinch of zinc dust and little	
solid NH4Cl. The mixture is cooled and	
filtered into Tollen's reagent.	

<u>5. Acetone</u> CH₃COCH₃

Physical characteristics and preliminary test :

1. State	: Liquid
2. Colour	: Colorless
3. Odour	: Pleasant ethereal odour
4. Miscibility	: Miscible with water
5. Litmus	: Neutral
6. Action of heat	: Completely Volatilises without leaving any residue
7. Ignition test	: Blue non-sooty flame

Special Tests

Experiment	Observation
1. Brady's test : A small amount of O.S. is	1. Yellow ppt. on mixing.
added to 2-4-dinitrophenylhydrazine solution	
2. Denige's test: An equal volume of Denige's	
reagent is added to an aqueous solution of O.S	2. Heavy white ppt. of double compound of
and the test tube is then kept in a boiling water	acetone and basic mercuric sulphate.
bath for few minutes.	
	3. A ruby red (reddish-purple) color develops

3. Legal's test: A few drops of a very dil. which disappears on warming but reappears solution of sodium nitroprusside is added to 2-on cooling.
3 ml of aqueous solution of O.S. followed by a few drops of dil. NaOH solution.
4. (a) Yellow ppt. of CHI₃ having characteristics sweeet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Yellow ppt. of CHI₃ having characteristics sweet smell.
(b) Solution instead of NaOH (Distinction from ethanol)



Physical characteristics and preliminary test :

1. State

2. Colour

3. Odour

4. Miscibility

- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Liquid
- : Colorless when freshly distilled but turns brown on exposure to light and air
- : Characteristic aromatic smell
- : Immiscible with water but soluble
 - in dil. HCI
 - : Neutral
 - : Completely volatilises without leaving any residue
 - : Yellow sooty flame
- Special Tests

Experiment	Observation
1. Bleaching powder test: A few drops of	1. Purple-violet colouration.
bleaching powder solution is added to a dilute	
solution of O S. in HCl.	
2. Potassium dichromate test: A drop of O.S.	2. Intense blue colour develops.
is added to 5 to 6 drops of conc. H ₂ SO ₄ taken in	r
a spot plate and the mixture is stirred well with	
glass rod. Then a drop of K ₂ Cr ₂ O ₇ solution is	2 Prilliant scorlat rad dya
added to it.	5. Brittant scarlet-fed uye.
3. Diazocoupling test: 5-6 drops of O.S are	
dissolved in dil. HCl in a test tube and cooled	
in ice-water. Then 3-4 drops of very dilute ice	
cold solution of sodium nitrite are added to it.	
Then the solution is added dropwise to 2 mL of	

ice-cold alkaline solution of β -naphthol; finally solution should remain mild alkaline	
---	--



Physical characteristics and preliminary test :

1.	State	
2.	Colour	

- Odour
 Miscibility
- 5. Litmus
 6. Action of heat
- 7. Ignition test

- : Liquid
- : Colorless when freshly distilled but turns brown on exposure to light and air
- : Characteristic bad smell
- : Immiscible with water but soluble in dil.HCI
- : Neutral
- : Completely volatilises without
- leaving any residue
- : Yellow sooty flame

Special Tests

Experiment	Observation
1. Bleaching powder test: A few drops of	1. Purple-violet colouration.
bleaching powder solution is added to a dilute	
solution of O S. in HCl.	
2. Potassium dichromate test: A drop of O.S.	2. Intense blue colour develops.
is added to 5 to 6 drops of conc. H ₂ SO ₄ taken in	I I I I I I I I I I I I I I I I I I I
a spot plate and the mixture is stirred well with	
glass rod. Then a drop of K ₂ Cr ₂ O ₇ solution is	2 A colour appears at this stage
added to it.	S. A colour appears at this stage.
3. Nitrous acid test: 5-6 drops of O.S are	(No red dye after addition to alkaline 2-
dissolved in dil. HCl in a test tube and cooled	naphthol solution)
in ice-water. Then 3-4 drops of very dilute ice	(distinction from aniline)
cold solution of sodium nitrite are added to it.	
4. Malachite green test: See test of	4. An intense green colouration
benzaldehyde below. Use benzaldehyde and	
this sample.	

8. Benzaldehyde C₆H₅CHO

Physical characteristics and preliminary test :

1. State

- 2. Colour 3. Odour
- 4. Miscibility
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

: Liquid

- : Colorless
- : Characteristic smell of bitter
- almonds : Immiscible with water
- : Neutral
- : Volatilises
- : Yellow sooty flame

Special Tests

Experiment	Observation
1. 2,4-D.N.P test: A few drops of 2,4-D.N.P	1. Reddish yellow ppt. forms immediately
solution is added to 1 mL alcoholic solution	simply on mixing.
given sample.	
2. Tollen's test: An alcoholic solution of the	2. Black ppt.
given sample is added to 2 mL of Tollen's	
reagent, warmed gently.	3. A intense green coloration
3. Malachite green test: 0.5 ml of O.S is	
heated with 1 ml of dimethyl aniline and a small	
bit of anhydrous ZnCl ₂ in a dry test tube for one	
minute. The leucobase produced is oxidised	
with lead dioxide in a solution of acetic acid	
and excess conc. HCl is added.	

Solid compounds

<u>1. Oxalic acid</u>

[(COOH)2]

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystaline
3. Odour	: Odourless
4. Solubility	: Soluble in cold water
5. Litmus	: Blue litmus turns red
6. Action of heat	: Volatilises completely without charring when heated slowly and

a sublimation is formed but decomposed with evolution of gases when heated strongly

7. Ignitio	n test
------------	--------

: Blue non-sooty flame

Experiment	Observation
1. H ₂ SO ₄ test: A small amount of O.S. is	1. Lime water turns milky.
heated with conc. H ₂ SO ₄ and the evolved	
gas is passed into lime water.	
2. Soda-lime test : A small amount of O.S.	2. The gas burns at the mouth of the test
is heated in a hard glass test tube with	tube.
sodalime.	
3. Calcium chloride test: Calcium	3. A white ppt. forms immediately.
chloride solution is added to neutral	
solution of 0.5. The ppt is divided into two	
(i) Dil acatic acid is added to one portion	(i) The ppt does not dissolve.
(1) Diff. acetic acid is added to one portion.	(Distinction from tartaric and succinic
	acids)
(ii) Dil. HCl or dil. HNO ₃ is added to	(ii) The ppt. dissolves
another portion.	4. A white crystalline ppt.
4. Silver nitrate test : AgNO ₃ solution is	
added to neutral solution of O.S. The ppt.	(i) The ppt. dissolves
is divided into three parts :	(ii) The ppt. dissolves
(il Dil. NaOH is added to one portion. (ii)	(iii) Metallic silver is formed
Dil. HNO ₃ is added to another portion.	()
(iii) The third portion is warmed, dried and	5 White ppt in the cold which does not
heated strongly in a dry test tube	dissolve
5. Denige's test (C.T.): Denige's reagent is	6 The pink colour of permanganate
added to neutral solution or aqueous	disappears.
solution of O.S.	
6. Potassium permanganate test: Alittle	
of dil. H ₂ SO ₄ is added to a solution of O.S or	
neutral soln. of O.S. The solution is warmed	
and dil. KMnO4 solution is added drop by	
drop with shaking.	

2. Succinic acid

Physical characteristics and preliminary test :

1. State: Solid2. Colour: Colorless3. Texture: Crystaline3. Odour: Odourless4. Solubility: Soluble in cold water5. Litmus: Blue litmus turns red

6. Action of heat

: Melts and then boils giving off extremely irritating vapours. A sublimate is forme d at the cooler part of the test tube : Blue non-sooty flame

7. Ignition test

Experiment	Observation
1. H ₂ SO ₄ test: Few crystals of O.S. is	1. O.S. dissolves without charring. Slight
warmed conc H_2SO_4 and then heated	chairring occurs on strong heating and the
strongly	solution turn brown with evolution of SO
subligiy.	solution turn brown whitevolution of 302.
	2. A white ppt. soluble in acetic acld.
2. Calcium chloride test: An equal volume	(Distiction from oxalic macid)
of Calcium chloride solution is added to	
neutral solution of O.S. The mixture is	
shaken vigorously inner side of the test	
tube is scratched with glass rod and finally	3. Buff-colored ppt. soluble in dil. HCl.
boiled	4. No ppt. in the cold but the color of
3 Silver nitrate test \cdot AgNO ₂ solution is	permanganate presists.
added to neutral solution of OS	
4. Denige's test (C.1.): Denige's reagent is	
added to neutral solution or aqueous	5. A beautiful reddish green fluorescence,
solution of O.S. and then a drop of 2%	intesifies to a beautiful green
KMnO _{4 solution} is added.	fluorescence on addition of alkali
5. Fluorescein test: A small quantity of	
the sample ix mixed with resorcinol and 2	
/ 3 drops of c. sulfuric acid and heated.	
Then poured into large volume of water.	

3. Resorcinol

1. State :	Solid
2. Colour :	Colorless
3. Texture :	Crystaline
3. Odour :	Odourless
4. Solubility :	Soluble in cold water
5. Litmus :	feebly acidic
6. Action of heat	: Melts and volatiles
7. Ignition test :	Yellow sooty flame
Experiment	Observation
1. H ₂ SO ₄ test: Few crystals of O.S. is warmed with	1. O.S dissolve without charring.
conc. H_2SO_4 .	
2. Ferric chloride test: A few drops of FeCl ₃	2. Bluish-violet color.
solution is added to aqueous solution of O.S.	
3. Bromine writer test: Bromine writer is	3. White ppt.
added to aq. solution of O.S. and shaken.	4. Black ppt of silver mirror.
4. Tollen's test: Aqueous solution of O.S is added	
o Tollen's reagent and warmed.	5. An intense reddish-green fluorescein,
	intesifies to a beautiful green

5. Fluorescein test: A small quantity of O.S. is fluorescence on addition of alkali.
mixed with succinic acid / phathalic acid and 26. A brilliant reddish colouration is drops conc. Sulfuric acid and heated. Then poured developed in aq. layer.
into a large volume of water.
6. O.S. is mixed with 3 - 4 drops CHCl₃ and gently heated on a basin with NaOH solution.

	4. Urea	
Physical characteristics an	nd preliminary test :	
1. State	: Solid	
2. Colour	: Colorless	
3. Texture	: Crystaline	
3. Odour	: Odourless	
4. Solubility	: Soluble in cold water	
5. Litmus	: Neutral	
6. Action of heat	: Melts and gives off characteristic smell of NH ₃	
7. Ignition test	: Blue non-sooty flame	

Experiment	Observation
1. Nitrous acid test: Dil. HCl and dil. NaNO2	1. Effervescence with evolution of gas.
solution is added to aq. Solution of O.s in a test tube	
and shaken.	2. A white crystalline ppt.
2. Nitric acid test : 2 mL of conc. HNO ₃ is added to aqueous concentrated solution of O.S.	3. A white crystalline ppt.
3. Oxalic acid test: A concentrated solution of oxalic acid is added to conc. solution of O.S.	4. A pink or violet colour develops.
4. Biuret test: About 0 5 g of O.S. is heated gently [#] to melt. The heating is continued till the molten mass solidifies at once when test tube is taken out of	[#strong and rapid heaing leads to formation of cyanuric acid and then the test will fail]
flame. If there is no solidification the process is repeated. The residue is cooled and dissolved in I ml	
of NaOH solution by warming and shaking. The solution is cooled again and a drop or two of a very dilute copper sulphate solution is added to it.	

5. Benzoic acid

1. State	: Solid
2. Colour	: White
3. Texture	: Plate-shaped or needle shaped Crystal
3. Odour	: Odourless
4. Solubility	: Insoluble in cold water but Soluble in hot water
5. Litmus	: Blue litmus turns red
6. Action of heat	: Melts and sublimates
7. Ignition test	: Yellow sooty flame

Experiment	Observation
1. H ₂ SO ₄ test: A small amount of O.S. is warmed with conc. H_2SO_4 .	1. Dissolves without charring.
 2. Soda-lime test: A small amount of O.S. is heated in a hard glass test tube with sodalime. 3. Denige's test: Denige's reagent is added to neutral solution of O.S. 	 Characteristics smell of benzene. White ppt. dissolves on boiling but appears on cooling. (a) a buff colored ppt.
 4. Ferric chloride test: (a) A drop or two of freshly prepared solution of Fecl₃ is added to neutral solution of O.S. (b) Dilute HCl is added to it. 	 (a) a bull-colored ppt. (Distinction from salicylic acid) (b) buff-colored ppt dissolves with appearance of white ppt.
5. Esterification test: About 0.1 g of O.S and 2 mL of dehydrated alcohol are taken in a clean test tube and shaken well. Then a few drops of conc. H_2SO_4 is added and the test tube is waremd in a hot water bath. The test tube is cooled and content of the test tube is poured into 50 mL of dilute sodium cxarbonate solution taken in a beaker.	5. Characteristics fruity smell of ester.

6. Salicylic acid

- 2. Colour
- 3. Texture
- Odour
 Solubility
- .
- 5. Litmus
- 6. Action of heat
- 7. Ignition test

- : Solid
- : White
- : needle shaped Crystal
- : Odourless
- : Sparingly soluble in cold water but readily Soluble in hot water
- : Blue litmus turns red
- : Melts and sublimes
- : Yellow sooty flame

Experiment	Observation
 H₂SO₄ test: A small amount of O.S. is warmed with conc. H₂SO₄. Soda-lime test: A small amount of O.S. is boated in a boat along test tube with 	 The solid dissolved, charring occurs after some time, solution darkens and gas evolve. Characteristics smell of phenol.
 as heated in a hard glass test tube with sodalime. 3. Denige's test: Denige's reagent is added to neutral solution of O.S. 4. Ferric chloride test: A drop or two of freshly prepared solution of Fecl₃ is added to neutral solution of O.S. (a) Dilute HCl is added to it. (b) acetic acid is added to another part. 5. Oil of wintergreen test: About 0.1 g of O.S and 1 mL of methanol and few drops of conc. H₂SO₄ are warmed in aclean dry test tube. 	 3. White ppt. dissolves on boiling but appears on cooling. 4. An intense violet color (<i>Distinction from benzoic acid</i>) (a) The color discharges (b) The color persists 5. A characteristics pungent fragnant odor which intensifies on pouring the mixture to dil. Sodium carbonate solution.

7. Tartaric acid

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystalline
3. Odour	: Odourless
4. Solubility	: Soluble in cold water
5. Litmus	: Blue litmus turns red
6. Action of heat	: Chars and gives off smell of burnt sugar
7. Ignition test	: Blue non-sooty flame

Experiment	Observation
1. H ₂ SO ₄ test: A small amount of O.S. is	1. charrs immediately.
warmed with conc. H_2SO_4 .	
2. Calcium chloride test: Excess of	2. A white crystalline ppt. soluble in hot
calcium chloride solution is added to	dilute acetic acid.
neutral solution of O.S. shaken and inside	
of the test tube is scratched with a glass rod.	
4. Silver nitrate test : (a) AgNO ₃ solution	4 (a) A white crystalline ppt (b)
is added to neutral solution of O.S. (b)	shining silver mirror forms at the inner
Addition of AgNO3 solution is continued	of the test tube.
till the precipitation is complete. Then dil.	

NH4OH is added dropwise with constant	(Distinction from oxalic, citric, succinic
shaking till the ppt. almost and not	acids)
completely dissolved as revealed by slight	
turbidity.	
5. Denige's test (C.T.): Denige's reagent is	5 Permanganate color discharges
added to neutral solution or aqueous	immediately without turbidity
solution of O.S.	minediatery without turblenty.

8. Glucose

Physical characteristics and preliminary test :

1. State	: Solid
2. Colour	: Colorless
3. Texture	: Crystalline
3. Odour	: Odourless
4. Solubility	: Soluble in cold water
5. Litmus	: Neutral
6. Action of heat	: Charrs and gives off smell
	of burnt sugar
7. Ignition test	: Blue non-sooty flame
Experiment	Observation

H₂SO₄ test: A small amount of O.S. is 1. No charring in cold but darkens on warmed with conc. H₂SO₄.
 Sodium hydroxide test: 2 ml of 3%

solution of NaOH is added to conc. solution². The solution turns first yellow and then of O.S. and the mixture is heated. Theraddish brown.

mixture is then acidified with dil. HNO₃.

3. Tollen's test: 5 ml of aqueous solution of

O.S. is added to equal volume of Tollen's 3. Black or grey precipitate or shining reagent taken in a clean test tube and then the silver mirror.

test tube is placed in a beaker of boiling

water bath for a few minutes.

4. Fehling's test: The sample is added to a mixture of Fehling A and Felling B and ^{4.} Brick red precipitate.

heated for 1 minute.

5. Lead acetate test: Lead acetate solution is added to solution of o.s and the mixture is ^{5.} The white ppt. turns salmon-pink in color.

boiled for few seconds. Then dil. NH₄OH

solution is added dropwise till just sufficient to produce a permanat white ppt. then the mixture is boiled again.

SEC1P: Pharmaceutical Chemistry

1. Preparation of Aspirin and its analysis Principle:

Acetylsalicylic acid or aspirin is acetyl derivative of salicylic acid. It is prepared by acetylation reaction of salicylic acid using acetyl chloride or acetic anhydride under acidic condition.

Reaction:



Chemicals required:

Salicylic acid : 5 gm Acetyl Chloride : 9 ml Conc. H₂SO₄: 5 ml Ice cold water: 100 ml Ethanol: 10 mL FeCl₃ solution

Procedure:

1. 5 gm of salicylic acid was taken in a 200 ml of conical flask and to this add 9 ml of acetyl chloride was added and they were mixed properly with the help of a glass rod.

2. To this, a few drops of conc. H_2SO_4 was added and mixing was done with the help of a glass rod.

3. After that 50 mL of ice cold water was added to the reaction mixture, scratch the sides of flask with glass rod and filter the solution.

4. The residue was dried in hot oven and, weighed. and recrystallized from ethanol.

5. 10 mg of crude aspirin was taken in a test tube and dissolve it in about 1 ml of alcohol. 1 drop of 1% FeCl₃ solution was added to it. Absence of blue or green colour indicate the absence of unreacted salicylic acid.

6. Recrystalization was done by taking half of the crude sample and dissolve it in minimum quantity of ethanol and to this few drops warm water was added until a turbidity persist. The solution was then cooled, filtered and the filtrate was kept for crystallization.

7. The crystals were collected by filtration, washed and dried. Then weight was measured.

8. Melting point of the crystalized sample was checked.

Result:

Weight of crude product:

Weight of the crystalized product:

..... g of acetylsalicylic acid was obtained from 5 g of salicylic acid.

Theoretical yield = 6.52 g

% Yield = (Yield obtained \times 100)/Theoretical Yeild

. The melting point of aspirin was found to be = $\dots^{\circ}C$.

2. Preparation of magnesium bisilicate (Antacid):

Magnesium silicate is manufactured by the precipitation reaction between sodium silicate and a watersoluble magnesium salt such as magnesium chloride, magnesium nitrate or magnesium sulfate. Magnesium silicate powders were prepared by a hydrothermal method. Na₂SiO₃ with the volume of 0.9 mL was dissolved in water, while 1.48 g of Mg(NO₃)₂.6H₂O was dissolved in the mixture of propylene glycol–400 and ethanol with the volume ratio of 3 to 1. After the two solutions were mixed, a white precipitate was formed. The aqueous suspension of the precipitate is filtered and the collected solid washed, dried. After the magnesium silicate has been filtered and dried, the prepared cake discharged from the dryer may then be used in tablets, granulations and the like or the magnesium silicate may be added directly as a constituent of a liquid antacid preparation. The composition of the precipitate depends on the ratio of the components in the reaction medium, the addition of the correcting substances, and the way in which they are precipitated. The molecular formula is typically written as MgO: XSiO₂, where X denotes the average mole ratio of SiO₂ to MgO. The product is hydrated and the formula is sometimes written MgO: XSiO₂•H₂O to show the water of hydration. Although magnesium silicate is of variable composition, the molar ratio of MgO to SiO₂ is approximately 2:5.

The reaction involved in the described process may be typically illustrated by the following equation:

Wherein x and y depend upon the type silica fed into the processing system.

